# GGT BR

IFCC


*Enzymatic colorimetric method*

KINETIC

<b>REF CT10192</b> 3 x 50 mL	<b>REF CT10195</b> 12 x 50 mL
<b>CONTENTS</b> R1.Reagent 3 x 40 mL R2.Reagent 3 x 10 mL	<b>CONTENTS</b> R1.Reagent 12 x 40 mL R2.Reagent 12 x 10 mL
For <i>in vitro</i> diagnostic use only	

## SUMMARY

Gamma-glutamyltransferase ( $\gamma$ -GT) catalyzes the transfer of a  $\gamma$ -glutamyl group from  $\gamma$ -glutamyl-3-carboxy-4-nitroanilide to glycylglycine with the formation of L- $\gamma$ -glutamyl-glycylglycine and 5-amino-2-nitrobenzoate. The amount of 5-amino-2-nitrobenzoate formed, monitored kinetically at 405 nm, is proportional to the enzyme activity present in the sample.<sup>1</sup>

## REAGENTS

**R1. Buffer.** TRIS 133 mmol/L, glycylglycine 138 mmol/L.

**R2. Substrate.** L- $\gamma$ -Glutamyl-3-carboxy-4-nitroanilide 23 mmol/L.

## PREPARATION

The Reagents are ready-to-use.

## STORAGE AND STABILITY

Store at 2-8°C. The Reagents are stable until the expiry date stated on the label. After daily use stored tightly closed and protected from light. On board the reagents are stable 30 days. Discard the reagent if the blank presents an absorbance over 1.200 at 405 nm against distilled water.

## SAMPLE COLLECTION

Serum or EDTA plasma free of hemolysis. Fluoride, citrate and oxalate inhibit  $\gamma$ -GT activity<sup>2</sup>.

The enzyme in the sample is stable for at least 1 week at 2-8°C and for at least 2 months when frozen.

## INTERFERENCES<sup>6</sup>

- Lipemia (intralipid >2.5 g/L) may affect the results.
- Bilirubin (> 10 mg/dL) may affect the results.
- Hemoglobin (> 8 g/L) may affect the results.
- Other drugs and substances may interfere<sup>6</sup>.

## INSTRUMENTATION AND MATERIALS

- LIDA analyzer.
- Laboratory equipment.
- iso-Clean Solution. Ref. CT18002.
- Multicalibrator CC/H 10x3 mL Ref. CT19750

## AUTOMATED PROCEDURE

A graphic display pictures the specific sets corresponding to the technical application outlined for this test.

Any new application, to the instrument should be validated to confirm that results meet the analytical performance of the method.

It is recommended to validate periodically the instrument.

## CALIBRATION

No calibration is required. Activity is calculated against theoretical factor.

## CALCULATION

Patient values are calculated automatically by theoretical factor.

## RESULTS

Samples with  $\Delta A/\text{min}$  exceeding 0.200 at 405 nm should be diluted 1:10 with saline and assayed again. Multiply the results by 10.

If results are to be expressed as SI units apply:  $U/L \times 16.67 = \text{nkat/L}$

## EXPECTED VALUES<sup>4</sup>

Serum, plasma

Temperature	37°C	30°C	25°C
Men	10-50 U/L (167-834 nKat/L)	7-35 U/L (117-538 nKat/L)	5-25 U/L (83-417 nKat/L)
Women	8-35 U/L (133-583 nKat/L)	6-25 U/L (100-417 nKat/L)	6-25 U/L (83-300 nKat/L)

It is recommended that each laboratory establishes its own reference range.

## QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) with assayed values handled as unknowns. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

## DIAGNOSTIC CHARACTERISTICS

Gamma-glutamyl transferase is the most sensitive enzymatic indicator of hepatobiliary disease available. It is highest in cases of *intrahepatic* or *posthepatic biliary obstruction*, and is more sensitive than alkaline phosphatase in detecting *obstructive jaundice*, *colangitis*, and *cholecystitis*.

Elevated levels are noted in the sera of patients with *alcoholic cirrhosis* and from people who are heavy drinkers. The enzyme levels are important in detecting alcohol induced liver disease correlating well with the duration of the drug action.

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

## PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

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