# GLUCOSE MR


Enzymatic colorimetric method

ENDPOINT

<b>REF CT10200</b> 2 x 40 mL	<b>REF CT10202</b> 6 x 40 mL	<b>REF CT10205</b> 10 x 60 mL
<b>CONTENTS</b> R1.Reagent 2 x 40 mL	<b>CONTENTS</b> R1.Reagent 6 x 40 mL	<b>CONTENTS</b> R1.Reagent 10 x 60 mL
For <i>in vitro</i> diagnostic use only		

## SUMMARY

In the Trinder reaction<sup>1,2</sup>, the glucose is oxidized to D-gluconate by the glucose oxidase (GOD) with the formation of hydrogen peroxide. In the presence of peroxidase (POD), a mixture of phenol and 4-aminoantipyrine (4-AA) is oxidized by hydrogen peroxide, to form a red quinoneimine dye proportional to the concentration of glucose in the sample.

## REAGENT

**R1.** Phosphate buffer 100 mmol/L pH 7.5, glucose oxidase > 10 KU/L, peroxidase > 2 KU/L, 4-aminoantipyrine 0.5 mmol/L, phenol 5 mmol/L.

## PREPARATION

The Monoreagent is ready-to-use.

## STORAGE AND STABILITY

Store at 2-8°C. The Reagent is stable until the expiry date stated on the label. After daily use stored tightly closed and protected from light.

Discard the reagent if presents an absorbance over 0.100 at 500 nm against distilled water.

On board the reagent is stable 30 days.

## SAMPLE COLLECTION

Serum or heparin plasma free of hemolysis.

Glucose is stable up to 24 hours at 2-8°C when serum or plasma is separated within 30 minutes after collection.

## INTERFERENCES

- Lipemia (intralipid) may affect the results.
- Bilirubin (> 10 mg/dL), hemoglobin (> 1 g/L), may affect the results.
- Other drugs and substances may interfere<sup>4</sup>.

## INSTRUMENTATION AND MATERIALS

- LIDA analyzer.
- Laboratory equipment.
- iso-Clean Solution. Ref. CT18002.
- Multicalibrator CC/H 10x3 mL Ref. CT19750

## AUTOMATED PROCEDURE

A graphic display pictures the specific sets corresponding to the technical application outlined for this test.

Any new application, to the instrument should be validated to confirm that results meet the analytical performance of the method.

It is recommended to validate periodically the instrument.

## CALIBRATION

Recalibrate weekly, when a new lot of reagent is used, when control recovery falls out of the expected range or when adjustments are made to the instrument. Two point calibration is recommended (S1: NaCl 9 g/L and S2: Calibrator). A reagent blank should be run daily before sample analysis.

## RESULTS

Samples with concentrations higher than 500 mg/dL should be diluted 1:4 with saline and assayed again. Multiply the results by 4.

If results are to be expressed as SI units apply: g/dL x 0.0555 = mmol/L

## EXPECTED VALUES<sup>5</sup>

Serum, plasma (fasting)

Adults	70-105 mg/dL (3.89-5.83 mmol/L)
Children	60-110 mg/dL (3.33-6.11 mmol/L)
Newborns	40-60 mg/dL (2.22-3.33 mmol/L)

It is recommended that each laboratory establishes its own reference range.

## QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) with assayed values handled as unknowns.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

## DIAGNOSTIC CHARACTERISTICS

The glucose level in the blood is maintained by diet uptake and regulatory hormones such as insulin, glucagon, or epinephrine.

An abnormal increase in blood glucose level, referred as *hyperglycemia*, can be associated with diabetes mellitus and hyperactivity of thyroid, pituitary or adrenal glands.

An abnormal decrease beyond the fasting level, referred as *hypoglycemia*, is observed in cases of insulin overdose, insulin secreting tumors, mixedema, hypopituitarism, Addison's disease and conditions interfering with glucose absorption.

Glucose measurement in the blood is a key test to evaluate and diagnose any carbohydrate-related disorder.

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

## PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

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