IgG at

Turbidimetric method

REF CT31800 2 x 40 mL	REF CT31802 6 x 40 mL
CONTENTS R1.Reagent 2 x 40 mL	CONTENTS R1.Reagent 6 x 40 mL
For <i>in vitro</i> diagnostic use only	

SUMMARY

IgG at is a quantitative turbidimetric assay^{1,2} for the measurement of IgG in human serum or plasma.

Anti-human IgG antibodies form insoluble complexes when mixed with samples containing IgG. The scattering light of the immunocomplexes depends of the IgG concentration in the patient sample, and can be quantified by comparison from a calibrator of known IgG concentration.

REAGENTS

R1 IgG at. Goat antibodies anti-human IgG, tris buffer 20 mmol/L, pH 8.2.

Precautions: The reagent contains sodium azide 0.95 g/L. Avoid any contact with skin or mucous.

PREPARATION

R1. Ready to use.

Calibration curve:

Prepare two fold dilutions of the Calibrator (1/1, 1/2, 1/4, 1/8) using NaCl 9 g/L as diluent. Use NaCl 9 g/L as Calibrator 0. Calculate the concentration of each IgG Calibrator dividing the Calibrator concentration by the corresponding dilution factor.

STORAGE AND STABILITY

1.  The reagents will remain stable until the expiration date printed on the label. For optimal stability store **tightly closed** at 2-8°C. Do not use the reagents after the expiration date.
2. Presence of particles, turbidity and/or the absorbance of blank reagent > 0.3 at 540 nm are sign of deterioration.

SAMPLE COLLECTION

Fresh serum and EDTA or heparinized plasma.

IgG in serum or plasma is stable 7 days at 2-8°C or 3 months at -20°C.

Samples with presence of fibrin should be centrifuged before testing. Highly hemolyzed or lipemic samples are not suitable for testing.

INTERFERENCES

- Bilirubin (40 mg/dL) and rheumatoid factors (300 IU/mL) do not interfere. Hemoglobin (8 g/L), and lipemia (10 g/L) may affect the results. Other substances may interfere⁵.

ADDITIONAL EQUIPMENT

- Analyzer LIDA.

AUTHOMATIC TECHNIQUE

For automatic assays the application method is included.

Any application to an instrument should be validated to demonstrate that results meet the performance characteristics of the method.

It is recommended to validate periodically the instrument.

CALIBRATION

Recalibrate weekly, when a new lot of reagent is used, when control recovery falls out of the expected range or when adjustments are made to the instrument. A reagent blank should be run daily before sample analysis.

REFERENCE VALUES

Adults³: 700 – 1600 mg/dL

Newborn⁴: 299 – 852 mg/dL

It is recommended that each laboratory establishes its own reference range.

QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) with assayed values handled as unknown. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

DIAGNOSTIC CHARACTERISTICS³⁻⁵

Quantification of immunoglobulins in serum is important for the diagnosis of primary or secondary immunodeficiency, for monitoring of immunoglobulin therapy, and for following the clinical course of multiple myeloma⁴.

IgG is the most important immunoglobulin produced by plasma cells, and represents about 75% of the total immunoglobulins⁴.

Congenital and acquired immunodeficiency are the most important causes of IgG deficit^{3,4}.

Polyclonal hyperimmunoglobulinemia (normal response to infections), hepatitis, cirrhosis as well as autoimmune diseases are other causes that increase IgG concentration.

Increments of monoclonal IgG (paraproteins) are found in multiple myeloma, lymphocytic leukemia, and Waldenström macroglobulinemia⁴.

PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

NOTES

1. The linearity limit depends on the sample/reagent ratio, as well as the analyzer used. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
2. Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

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