Tranferrin at

Turbidimetric method

REF CT32100 2 x 40 mL CONTENTS R1.Reagent 2 x 40 mL	REF CT32105 6 x 40 mL CONTENTS R1.Reagent 6 x 40 mL
For <i>in vitro</i> diagnostic use only	

SUMMARY

Transferrin at is a quantitative turbidimetric assay^{1,2} for the measurement of transferrin (TRF) in human serum or plasma.

Anti-human TRF antibodies form insoluble complexes when mixed with samples containing TRF. The scattering light of the immunocomplexes depends of the TRF concentration in the patient sample, and can be quantified by comparison from a calibrator of known TRF concentration.

REAGENTS

R1 Transferrin at. Goat antibodies anti-human TRF, tris buffer 20 mmol/L, pH 8.2.

Precautions: The reagent contains sodium azide 0.95 g/L. Avoid any contact with skin or mucous.

PREPARATION

R1. Ready to use.

Calibration curve:

Prepare two fold dilutions of the Calibrator (1/1, 1/2, 1/4, 1/8) using NaCl 9 g/L as diluent. Use NaCl 9 g/L as Calibrator 0. Calculate the concentration of each Transferrin Calibrator dividing the Calibrator concentration by the corresponding dilution factor.

STORAGE AND STABILITY

1.  The reagents will remain stable until the expiration date printed on the label. For optimal stability store **tightly closed** at 2-8°C. Do not use the reagents after the expiration date.
2. Presence of particles, turbidity and/or the absorbance of blank reagent > 0.3 at 540 nm are sign of deterioration.

SAMPLE COLLECTION

Fresh serum and EDTA or heparinized plasma.

TRF in serum or plasma is stable 7 days at 2-8°C or 3 months at -20°C.

Samples with presence of fibrin should be centrifuged before testing. Highly hemolyzed or lipemic samples are not suitable for testing.

INTERFERENCES

- Bilirubin (40 mg/dL), hemoglobin (8 g/L), and rheumatoid factors (300 IU/mL) do not interfere. Lipemia (2.5 g/L) may affect the results. Other substances may interfere⁵.

ADDITIONAL EQUIPMENT

- Analyzer LIDA.

AUTHOMATIC TECHNIQUE

For automatic assays the application method is included.

Any application to an instrument should be validated to demonstrate that results meet the performance characteristics of the method.

It is recommended to validate periodically the instrument.

CALIBRATION

Recalibrate weekly, when a new lot of reagent is used, when control recovery falls out of the expected range or when adjustments are made to the instrument.

A reagent blank should be run daily before sample analysis.

REFERENCE VALUES

Adults³: 200 – 360 mg/dL

Newborn⁴: 117 – 250 mg/dL

It is recommended that each laboratory establishes its own reference range.

QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) with assayed values handled as unknown. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

DIAGNOSTIC CHARACTERISTICS³⁻⁵

Transferrin is the principal plasma protein for the transport of iron. It is synthesized in the liver and transfers iron through the serum. The TRF molecule specifically binds to the Fe³⁺ forming the TRF-Fe³⁺, that goes through the plasma and carries iron to storage sites in the body.

Evaluation of plasma TRF levels is useful for the differential diagnosis of anemia, for monitoring its treatment and for assessing the nutritional status of a patient.

TRF level rises in the hypochromic anemia (iron deficiency), in pregnancy and during estrogen administration.

Low levels of TRF occur in inflammation and malignancy.

PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

NOTES

1. The linearity limit depends on the sample/reagent ratio, as well as the analyzer used. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
2. Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

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