# C4 at

## *Turbidimetric method*

<b>REF CT31500</b> 2 x 40 mL	<b>REF CT31502</b> 6 x 40 mL
<b>CONTENTS</b> R1.Reagent 2 x 40 mL	<b>CONTENTS</b> R1.Reagent 6 x 40 mL
For <i>in vitro</i> diagnostic use only	

### SUMMARY

C4 at is a quantitative turbidimetric assay<sup>1,2</sup> for the measurement of (C4) in human serum or plasma.

Anti-human C4 antibodies form insoluble complexes when mixed with samples containing C4. The scattering light of the immunocomplexes depends of the C4 concentration in the patient sample, and can be quantified by comparison from a calibrator of known C4 concentration.

### REAGENTS

**R1 C4 at.** Goat antibodies anti-human C4, tris buffer 20 mmol/L, pH 8.2.

**Precautions:** The reagent contains sodium azide 0.95 g/L. Avoid any contact with skin or mucous.

### PREPARATION

**R1.** Ready to use.

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Prepare two fold dilutions of the Calibrator (1/1, 1/2, 1/4, 1/8) using NaCl 9 g/L as diluent. Use NaCl 9 g/L as Calibrator 0. Calculate the concentration of each C4 Calibrator dividing the Calibrator concentration by the corresponding dilution factor.

### STORAGE AND STABILITY

1.  The reagents will remain stable until the expiration date printed on the label. For optimal stability store **tightly closed** at 2-8°C. Do not use the reagents after the expiration date.
2. Presence of particles, turbidity and/or the absorbance of blank reagent > 0.3 at 340 nm are sign of deterioration.

### SAMPLE COLLECTION

Fresh serum and EDTA or heparinized plasma.

C4 in serum or plasma is stable 7 days at 2-8°C or 3 months at -20°C. Samples with presence of fibrin should be centrifuged before testing. Highly hemolyzed or lipemic samples are not suitable for testing.

### INTERFERENCES

- Bilirubin (10 mg/dL) and rheumatoid factors (400 IU/mL) do not interfere. Hemoglobin (4 g/L) and lipemia (2.5 g/L) may affect the results. Other substances may interfere<sup>7</sup>.

### ADDITIONAL EQUIPMENT

- Analyzer LIDA.

### AUTHOMATIC TECHNIQUE

For automatic assays the application method is included.

Any application to an instrument should be validated to demonstrate that results meet the performance characteristics of the method.

It is recommended to validate periodically the instrument.

### CALIBRATION

Recalibrate weekly, when a new lot of reagent is used, when control recovery falls out of the expected range or when adjustments are made to the instrument. A reagent blank should be run daily before sample analysis.

### REFERENCE VALUES

Adults<sup>3</sup>: 10 – 40 mg/dL

Newborn<sup>4</sup>: 13 – 38 mg/dL

It is recommended that each laboratory establishes its own reference range.

### QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) with assayed values handled as unknown. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

### DIAGNOSTIC CHARACTERISTICS<sup>3-5</sup>

C4 is the complement component essential for classical pathway activation. Most individuals with C4 deficiency do not have problems with infection, suggesting that the alternative pathway can compensate for the lack in the classical pathway activation in removal of bacterial agents. Hepatic cells synthesize C4, although in less proportion may be synthesized by monocytes and other tissues.

Increased and decreased levels of C4 both have clinical significance.

Increased levels are closely related with acute-phase response (trauma, inflammatory process).

Decreased levels are related with genetic deficiency (autoimmune or collagen vascular disease, particularly Systemic Lupus Erythematosus), or acquired deficiency as a consequence of the consumption in immunocomplexes formation, autoimmune hemolytic anemia and autoimmune nephritis.

### PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

### NOTES

1. The linearity limit depends on the sample/reagent ratio, as well as the analyzer used. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
2. Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

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