PHOSPHORUS UV

INORGANIC


Ultraviolet method

ENDPOINT

REF CT10320	REF CT10322	REF CT10325
2 x 40 mL	6 x 40 mL	10 x 60 mL
CONTENTS	CONTENTS	CONTENTS
R1.Reagent 2 x 40 mL	R1.Reagent 6 x 40 mL	R1.Reagent 10 x 60 mL
For <i>in vitro</i> diagnostic use only		

SUMMARY

Inorganic phosphate reacts with ammonium molybdate in the presence of sulfuric acid to form a phosphomolybdic complex which is measured at 340 nm.^{1,2}

REAGENT

R1. Molybdate Reagent. Ammonium molybdate 0.40 mmol/L, sulfuric acid 210 mmol/L. **X₁ R:36/37/38 S:26/37/39.**

PREPARATION

The Reagent is ready-to-use.

STORAGE AND STABILITY

Store at 2-8°C. The Reagent is stable until the expiry date stated on the label. With time, the reagent may develop a light blue color which does not interfere with the assay. On board this reagent is stable 30 days.

Discard the reagent if the blank presents an absorbance over 0.600 at 340 nm against distilled water.

SAMPLE COLLECTION

Serum or heparinized plasma separated from cells as soon as possible, and urine (see Notes). Stable for 7 days at 2-8°C. Freeze for longer storage. Phosphorus in acidified samples of urine is stable for 6 months at 2-8°C.

INTERFERENCES

- Lipemia (intralipid) may affect the results.
- Bilirubin (< 10 mg/dL) does not interfere.
- Hemoglobin may affect the results.
- Other drugs and substances may affect the phosphorus values.³

INSTRUMENTATION AND MATERIALS

- LIDA analyzer.
- Laboratory equipment.
- iso-Clean Solution. Ref. CT18002.
- Multicalibrator CC/H 10x3 mL Ref. CT19750

AUTOMATED PROCEDURE

A graphic display pictures the specific sets corresponding to the technical application outlined for this test.

Any new application, to the instrument should be validated to confirm that results meet the analytical performance of the method.

It is recommended to validate periodically the instrument.

CALIBRATION

Recalibrate weekly, when a new lot of reagent is used, when control recovery falls out of the expected range or when adjustments are made to the instrument. Two point calibration is recommended (S1: NaCl 9 g/L and S2: Calibrator). A reagent blank should be run daily before sample analysis.

RESULTS

Samples with concentrations higher than 22 mg/dL should be diluted 1:2 with saline and assayed again. Multiply the results by 2.

If results are to be expressed as SI units apply: mg/dL x 0.323 = mmol/L

EXPECTED VALUES⁴

Serum, plasma

Children	4.0-7.0 mg/dL (1.29-2.26 mmol/L)
Men	2.5-4.5 mg/dL (0.81-1.45 mmol/L)
Women	1.50-6.8 mg/dL (0.48-2.19 mmol/L)

Urine

0.4-1.0 g/24-h (12.9-32.3 mmol/24-h)

It is recommended that each laboratory establishes its own reference range.

QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) with assayed values handled as unknowns. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

DIAGNOSTIC CHARACTERISTICS

Phosphorus and calcium metabolism are intertwined. In healthy persons, as serum calcium levels rise, those of phosphorus fall. Control of phosphorus levels is in part accomplished by regulation of renal excretion. However, fairly rapid fluctuations in serum inorganic phosphate can occur because the serum inorganic phosphate concentration is influenced by carbohydrate metabolism.

In *diabetes*, severe loss of phosphate is possible, since carbohydrate metabolism is deranged and phosphate tends to pass from the cell into extracellular fluid and then into plasma. It is then extracted and excreted by the kidney.

Increased levels are associated with *hypoparathyroidism*, during insulin treatment of *diabetic coma*, and with chronic *nephritis* rising as renal failure progresses.

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

NOTES

- Collect a 24-hor urine specimen into a plastic bottle containing 20 mL of 50% (v/v) HCl. Bring to 2 L with distilled water. Mix completely and test as described for serum.

PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

BIBLIOGRAPHY

1. Daly J.A. and Erthinghsausen G. Clin. Chem. 18 : 263 (1972).
2. Gamst, O. and Try, K. Scand. J. Clin. Lab. Invest. 40:486 (1980).
3. Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AACC Press, 2000.
4. Tietz. N.W. Clinical Guide to Laboratory Tests, 3rd Edition. W.B. Saunders Co. Philadelphia, PA. (1995).

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