MAGNESIUM MR

CALMAGITE


Colorimetric method

ENDPOINT

REF CT10300 2 x 40 mL CONTENTS R1.Reagent 2 x 40 mL
For <i>in vitro</i> diagnostic use only

SUMMARY

The method¹ is based on the specific binding of calmagite, a metallochromic indicator², and magnesium at alkaline pH with the resulting shift in the absorption wavelength of the complex. The intensity of the chromophore formed is proportional to the concentration of magnesium in the sample.

REAGENTS

R1. Chromogen. Calmagite 75 mmol/L, EGTA 60 mmol/L, amino-methyl-propanol 0.2 mol/L, KCl 0.2 mol/L, surfactant 0.05 % (w/v)

PREPARATION

The Reagent is ready to use.

STORAGE AND STABILITY

 Store at 2-25°C.

The reagent is stable until the expiry date stated on the label.

On board this reagent is stable 30 days. After daily use, stored **tightly closed and protected from light**.

Discard the mixture if the blank presents an absorbance above 0.600 at 520 nm against distilled water (see Notes).

SAMPLE COLLECTION

Serum or heparinized plasma free of hemolysis. Other anticoagulants (EDTA, oxalate and citrate) must not be used.

Magnesium in serum or plasma is stable for 10 days at 2-8°C. Freeze for longer storage.

INTERFERENCES³

- Lipemia (intralipid) may affect the results.
- Bilirubin (< 10 mg/dL) does not interfere.
- Hemoglobin (> 1 g/L) may affect the results.
- Other drugs and substances may interfere³.
- Many detergents and water supplies (see Notes).
- The interference by calcium is prevented by the use of EGTA into the buffer.

MATERIALS REQUIRED

- Analyzer LIDA.

AUTHOMATIC TECHNIQUE

For automatic assays the application is included.

Any application to an instrument should be validated to demonstrate that results meet the performance characteristics of the method.

It is recommended to validate periodically the instrument.

CALIBRATION

Recalibrate every day, when a new lot of reagent is used, when control recovery falls out of the expected range or when adjustments are made to the instrument. Two point calibration is recommended (S1: NaCl 9 g/L and S2: Calibrator). A reagent blank should be run daily before sample analysis.

RESULTS

Samples with concentrations higher than 10 mg/dL (4.12 mmol/L) should be diluted 1:2 with saline and assayed again. Multiply the results by 2.

If results are to be expressed as SI units apply: mg/dL x 0.412 = mmol/L

EXPECTED VALUES⁴

Serum, plasma

Children (2-12 years)	1.7-2.3 mg/dL (0.70-0.94 mmol/L)
Adults (12-60 years)	1.6-3.0 mg/dL (0.66-1.23 mmol/L)

It is recommended that each laboratory establishes its own reference range.

QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) with assayed values handled as unknowns. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

DIAGNOSTIC CHARACTERISTICS

Magnesium is considered an essential nutrient and a major intracellular cation. More than 50% of the total magnesium found in the body is complexed with calcium in the skeleton; however, only 1% of total body magnesium is found in the circulation. Between 60% and 70% of serum magnesium is free; the remaining percentage is bound to albumin, phosphate, citrate, and other ions.

Magnesium is also used as an *activator* (inorganic ion that is needed as a cofactor for an enzyme reaction) for more 300 enzymes, especially those involved in oxidative phosphorylation, glycolysis, cell replication, and protein synthesis.

The absorption of magnesium in the intestinal tract and its reabsorption in the renal tubules are directly associated with the body's need, balance and health.

Hypomagnesemia (an abnormal decrease in the level of serum magnesium) is usually associated with severe prolonged diarrhea, impairment of neuromuscular function, gastrointestinal malabsorption and alcoholism.

Hypermagnesemia (an abnormal elevation in the level of serum magnesium) is usually associated with renal glomerular failure, dehydration, severe diabetic acidosis, and Addison's disease.

PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

NOTES

- Most of the detergents and water softening products used in the labs contain chelating agents. A defective rinsing will invalidate the procedure. Keep the glassware acid washed and thoroughly rinsed at all times.
- Do not reuse old containers or mix fresh reagent with old reagent

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