# AMYLASE MR

Colorimetric enzymatic method


KINETIC

<b>REF CT10060</b> 2 x 40 mL	<b>REF CT10062</b> 6 x 40 mL
<b>CONTENTS</b> R1.Reagent 2 x 40 mL	<b>CONTENTS</b> R1.Reagent 6 x 40 mL
For <i>in vitro</i> diagnostic use only	

## SUMMARY

In this direct method<sup>1,2</sup>  $\alpha$ -amylase catalyzes the hydrolysis of 2-chloro-p-nitrophenyl- $\alpha$ -D-maltotrioside (CNP-G3) substrate at pH 6.0 forming 2-chloro-p-nitrophenol (CNP) and free glycosides.

The reaction is monitored kinetically at 405 nm by the rate of formation of the colored CNP produced, proportional to the activity of the  $\alpha$ -amylase in the sample.

CNP-G2 = 2- Chloro-nitrophenyl- $\alpha$ -D-maltoside

G3 = Maltotriose

G = Glucose

## REAGENT

**R1.** MES buffer 50 mmol/L pH 6.0, calcium acetate 5 mmol/L, sodium chloride 300 mmol/L, sodium thiocyanate 450 mmol/L, CNP-G3 2.25 mmol/L.

## PREPARATION

The Monoreagent is ready-to-use.

## STORAGE AND STABILITY

 Store at 2-8°C. The Reagent is stable until the expiry date stated on the label. After daily use stored tightly closed and protected from light. On board the reagent is stable 30 days.

Discard the reagent if the blank presents an absorbance over 0.250 at 405 nm against distilled water.

## SAMPLE COLLECTION

Serum, heparinized plasma and urine.

Serum and plasma  $\alpha$ -amylase is stable for 30 days at 2-8°C.

Random urine samples should be clear and precipitate free for testing.

Check the pH. Urines with a pH < 5 may reduce the enzyme stability.

Stable for 10 days at 2-8°C.

## INTERFERENCES

- Lipemia (intralipid 20 g/L) does not interfere.
- Bilirubin (40 mg/dL) does not interfere.
- Hemoglobin (16 g/L) does not interfere.
- Other drugs and substances may interfere<sup>3,5</sup>.

## INSTRUMENTATION AND MATERIALS

- LIDA analyzer.
- Laboratory equipment.
- iso-Clean Solution. Ref. CT18002.
- Multicalibrator CC/H 10x3 mL. Ref. CT19750

## AUTOMATED PROCEDURE

A graphic display pictures the specific sets corresponding to the technical application outlined for this test.

Any new application, to the instrument should be validated to confirm that results meet the analytical performance of the method.

It is recommended to validate periodically the instrument.

## CALCULATION

Patient values are calculated automatically by theoretical factor.

## CALIBRATION

No calibration is required. Activity is calculated against theoretical factor.

## RESULTS

Samples with  $\Delta A/\text{min}$  exceeding 0.500 at 405 nm should be diluted 1:10 with saline and assayed again. Multiply the results by 10.

If results are to be expressed as SI units apply: U/L x 0.01667 =  $\mu\text{kat/L}$

## EXPECTED VALUES<sup>4</sup>

Serum, plasma

< 86 U/L (1.43 $\mu\text{kat/L}$ )
------------------------------------

Urine

< 470 U/L (7.83 $\mu\text{kat/L}$ )
-------------------------------------

It is recommended that each laboratory establishes its own reference range.

## QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) with assayed values handled as unknowns.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

## DIAGNOSTIC CHARACTERISTICS

Amylase activity tests in serum and urine are mainly used in the diagnosis of diseases of the pancreas and in the investigation of pancreatic function. Amylase is found chiefly in the saliva and in pancreatic tissue. Normally, small amounts of amylase are present in the blood, but with various forms of pancreatic disturbance large amounts of amylase are secreted into the blood by the pancreas.

*Increased* levels are found associated with acute pancreatitis, pancreatic duct obstruction, intra-abdominal diseases, mumps and bacterial parotitis.

A significant amount of the serum amylase is excreted in the urine, and as a result elevation of serum activity is reflected in the rise of urinary amylase activity. Urine amylase appears to be more frequently elevated, reaches higher levels, and persists for longer periods.<sup>4</sup>

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

## PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

## BIBLIOGRAPHY

1. Winn-Deen, E.S., David, H., Sigler, G, and Chavez, R. Clin. Chem. 34 : 2005.
2. International Federation of Clinical Chemistry (IFCC). Clin. Chem. Lab. Med. 36 : 185.
3. Young DS. Effects of drugs on clinical laboratory tests, 4<sup>th</sup> ed. AACC Press, 1995.
4. Tietz. Textbook of Clinical Chemistr, 2<sup>nd</sup> Edition. Burtis CA, Ashwood ER. WB Saunders Co., 1994.
5. Friedman and Young. Effects of disease on clinical laboratory tests, 4<sup>th</sup> ed. AACC Press 2001.

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