CREATINE KINASE MB

CK-MB INHIBITED

Immunological UV method


KINETIC

REF CT10170 1 x 50 mL	REF CT10172 3 x 50 mL
CONTENTS R1.Reagent 1 x 40 mL R2.Reagent 1 x 10 mL	CONTENTS R1.Reagent 3 x 40 mL R2.Reagent 3 x 10 mL
For <i>in vitro</i> diagnostic use only	

SUMMARY

The major CK activity in normal serum is due to CK-MM and CK-MB isoenzymes, found mainly in skeletal and heart muscle. CK-BB is usually present in serum at low concentration.

Both, creatine kinase (CK) enzymes are dimers formed by the association of two sub-units from muscle (M) and nerve cells (B). The immunoinhibition from an specific antibody of both, MM sub-units and the single M sub-unit of CK-MB, allows the determination of the B sub-unit. The CK-B activity corresponding to half of CK-MB is measured by the increasing rate of absorbance resulting from the following coupled reactions.^{1,2}

REAGENTS

R1. Buffer/Glucose/NAC. Imidazol buffer 100 mmol/L pH 6.7, glucose 20 mmol/L, NAC 20 mmol/L, magnesium acetate 10 mmol/L, NADP 2.5 mmol/L, HK ≥ 4 KU/L, EDTA 2 mmol/L.

R2. Substrate/Coenzymes. CP 30 mmol/L, AMP 5 mmol/L, ADP 2 mmol/L, di(adenosine-5') pentaphosphate 10 μ mol/L, G6P-DH ≥ 1.5 KU/L. Sufficient CK-M human antibody to inhibit ≥ 3000 U/L of CK-MM at 37°C.

PREPARATION

Working reagent. Mix 4 mL of R1 + 1 mL of R2.

STORAGE AND STABILITY

Store at 2-8°C. The Reagents are stable until the expiry date stated on the label. After daily use stored tightly closed and protected from light. On board are stable 15 days. Discard the reagents if presents an absorbance over 0.400 at 340 nm against distilled water.

SAMPLE COLLECTION

Serum. Stable 8 days at 2-8°C or 1 month at -20°C. Chill the samples as rapidly as possible after collection.

Moderately or severely hemolyzed specimens are unsatisfactory for testing, as well as plasmas containing EDTA, heparin, citrate, or fluoride since they may produce unpredictable reaction rates.

INTERFERENCES

- Lipids (intralipid >1 g/L) may affect the results.
- Bilirubin (<20 mg/dL) does not interfere.
- Hemoglobin (>8 g/L) may affect the results.
- The presence of above normal concentrations of CK-BB or adenilatekinasa, and of macro or mitochondrial CK may affect the results⁶.
- Other drugs and substances may interfere⁸⁻⁹.

INSTRUMENTATION AND MATERIALS

- LIDA analyzer.
- Laboratory equipment.
- iso-Clean Solution. Ref. CT18002.
- Multicalibrator CC/H 10x3 mL Ref. CT19750

AUTOMATED PROCEDURE

A graphic display pictures the specific sets corresponding to the technical application outlined for this test. Any new application, to the instrument

CALIBRATION

No calibration is required. Activity is calculated against theoretical factor.

CALCULATION

Patient values are calculated automatically by theoretical factor.

RESULTS

Percentage of CK-MB activity:

$$\frac{\text{CK-MB Activity}}{\text{Total CK Activity}} \times 100 = \% \text{ CK-MB activity}$$

EXPECTED VALUES

The normal range of CK-MB activity in the serum of adults is considered to be 2.0 U/L to 19.5 U/L at 37°C. Newborns, infants, and children have higher serum CK-MB values than adults.

A ratio between CK-MB and total CK activities above 4% should be considered suspicious, and above 10% consistent with acute myocardial infarction.⁴ When the three conditions stated below are met the probability of myocardial injury is very high.⁵

Condition	Test temperature (37°)	
	1	CK Men
	CK Women	>167 U/L (2.78 μ kat/L)
2	CK-MB	>24 U/L (0.40 μ kat/L)
3	% CK-MB	6-25%

It is recommended that each laboratory establishes its own reference range.

QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) with assayed values handled as unknowns. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

DIAGNOSTIC CHARACTERISTICS

The major medical application of CK-MB assays is in adults for the diagnosis of acute myocardial infarction (AMI) or, the differentiation of myocardial injury from skeletal muscle damage. Creatine kinase-MB is considered one of the best laboratory indicators of AMI. CK-MB determination within the proper time frame after infarction is most critical being the useful interval ("window") from about 10 to 24 hour after the infarction. Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

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