# IgM at

*Turbidimetric method*

<b>REF CT31900</b> 2 x 40 mL	<b>REF CT31902</b> 6 x 40 mL
<b>CONTENTS</b> R1.Reagent 2 x 40 mL	<b>CONTENTS</b> R1.Reagent 6 x 40 mL
For <i>in vitro</i> diagnostic use only	

## SUMMARY

*IgM at* is a quantitative turbidimetric assay<sup>1,2</sup> for the measurement of IgM in human serum or plasma.

Anti-human IgM antibodies form insoluble complexes when mixed with samples containing IgM. The scattering light of the immunocomplexes depends of the IgM concentration in the patient sample, and can be quantified by comparison from a calibrator of known IgM concentration.

## REAGENTS

**R1 IgM at.** Goat antibodies anti-human IgM, tris buffer 20 mmol/L, pH 8.2.

**Precautions:** The reagent contains sodium azide 0.95 g/L. Avoid any contact with skin or mucous.

## PREPARATION

**R1.** Ready to use.

### Calibration curve:

Prepare two fold dilutions of the Calibrator (1/1, 1/2, 1/4, 1/8) using NaCl 9 g/L as diluent. Use NaCl 9 g/L as Calibrator 0. Calculate the concentration of each IgM Calibrator dividing the Calibrator concentration by the corresponding dilution factor.

## STORAGE AND STABILITY

1.  The reagents will remain stable until the expiration date printed on the label. For optimal stability store **tightly closed** at 2-8°C. Do not use the reagents after the expiration date.
2. Presence of particles, turbidity and/or the absorbance of blank reagent > 0.3 at 340 nm are sign of deterioration.

## SAMPLE COLLECTION

Fresh serum and EDTA or heparinized plasma.

IgM in serum or plasma is stable 7 days at 2-8°C or 3 months at -20°C.

Samples with presence of fibrin should be centrifuged before testing. Highly hemolyzed or lipemic samples are not suitable for testing.

## INTERFERENCES

- Bilirubin (20 mg/dL) and rheumatoid factors (400 IU/mL) do not interfere. Hemoglobin (4 g/L), and lipemia (1.25 g/L) may affect the results. Other substances may interfere<sup>5</sup>.

## ADDITIONAL EQUIPMENT

- Analyzer LIDA.

## AUTHOMATIC TECHNIQUE

For automatic assays the application method is included.

Any application to an instrument should be validated to demonstrate that results meet the performance characteristics of the method.

It is recommended to validate periodically the instrument.

## CALIBRATION

Recalibrate weekly, when a new lot of reagent is used, when control recovery falls out of the expected range or when adjustments are made to the instrument. A reagent blank should be run daily before sample analysis.

## REFERENCE VALUES

Adults<sup>3</sup>: 40 – 230 mg/dL

Newborn<sup>4</sup>: 5 – 30 mg/dL

It is recommended that each laboratory establishes its own reference range.

## QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) with assayed values handled as unknown. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

## DIAGNOSTIC CHARACTERISTICS<sup>3-5</sup>

Quantification of immunoglobulins in serum is important for the diagnosis of primary or secondary immunodeficiency, for monitoring of immunoglobulin therapy, and for following the clinical course of multiple myeloma<sup>4</sup>.

IgM is the only immunoglobulin synthesized by neonates. In adult serum, represents the 5 to 10% of the total immunoglobulins. It is a pentameric molecule, and its big size prevents its passage the extravascular spaces.

Congenital and acquired immunodeficiency are causes of IgM deficit<sup>3,6</sup>.

Polyclonal hyperimmunoglobulinemia (normal response to infections) increases IgM concentration, specially in primary viral infections and blood stream infections such as malaria.

IgM increases also are found in primary biliary cirrhosis and chronic active hepatitis.

Increments of monoclonal IgM (paraprotein) are found in proliferative disorders of plasma cells as Waldenström's macroglobulinemia<sup>4</sup>.

## PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

## NOTES

1. The linearity limit depends on the sample/reagent ratio, as well as the analyzer used. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
2. Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

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