URIC ACID MR


Enzymatic colorimetric method

ENDPOINT

REF CT10390 2 x 40 mL	REF CT10392 6 x 40 mL	REF CT10395 10 x 60 mL
CONTENTS R1.Reagent 2 x 40 mL	CONTENTS R1.Reagent 6 x 40 mL	CONTENTS R1.Reagent 10 x 60 mL
For <i>in vitro</i> diagnostic use only		

SUMMARY

Uric acid is oxidized by uricase to allantoin with the formation of hydrogen peroxide. In the presence of peroxidase (POD), a mixture of dichlorophenol sulphonate (DCBS) and 4-aminoantipyrine (4-AA) is oxidized by hydrogen peroxide to form a quinoneimine dye proportional to the concentration of uric acid in the sample.^{1,2}

REAGENT

R1. Phosphate buffer 100 mmol/L pH 7.8, uricase > 50 U/L, peroxidase > 1 KU/L, ascorbate oxidase > 0,1 KU/L, 4-aminoantipyrine 0.32 mmol/L, DCBS 2 mmol/L, non-ionic tensioactives 2 g/L (w/v). Biocides.

PREPARATION

The Monoreagent is ready-to-use.

STORAGE AND STABILITY

Store at 2-8°C. The Reagent is stable until the expiry date stated on the label. Discard the reagent if presents an absorbance over 0.100 at 520 nm against distilled water. After daily use stored tightly closed and protected from light. On board the reagent is stable 30 days.

SAMPLE COLLECTION

Whenever possible medication should be suspended 12 hours before sample collection. Hemolysis-free serum, EDTA or heparinized plasma (see Notes). Uric acid in serum or plasma is stable up to 5 days at 2-8°C and for 6 months at -20°C.

INTERFERENCES

- Lipemia (intralipid < 5 g/L) does not interfere.
- Bilirubin (< 10 mg/dL) does not interfere.
- Hemoglobin may affect the results.
- Other drugs and substances may interfere³.

INSTRUMENTATION AND MATERIALS

- LIDA analyzer.
- Laboratory equipment.
- iso-Clean Solution. Ref. CT18002.
- Multicalibrator CC/H 10x3 mL Ref. CT19750

AUTOMATED PROCEDURE

A graphic display pictures the specific sets corresponding to the technical application outlined for this test.

Any new application, to the instrument should be validated to confirm that results meet the analytical performance of the method.

It is recommended to validate periodically the instrument.

CALIBRATION

Recalibrate weekly, when a new lot of reagent is used, when control recovery falls out of the expected range or when adjustments are made to the instrument. A reagent blank should be run daily before sample analysis.

RESULTS

Samples with concentrations higher than 20 mg/dL should be diluted 1:5 with saline and assayed again. Multiply the results by 5.

If results are to be expressed as SI units apply: mg/dL x 59.5 = μmol/L

EXPECTED VALUES⁴

Serum, plasma

Men	3.5-7.2 mg/dL (208-428 μmol/L)
Women	2.6-6.0 mg/dL (155-357 μmol/L)

It is recommended that each laboratory establishes its own reference range.

QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) with assayed values handled as unknowns.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

DIAGNOSTIC CHARACTERISTICS

An abnormal increase in the level of uric acid in the circulation above 7.0 mg/dL (0.42 mmol/L) is referred to as *hyperuricemia*, being the gout the major form of the ailment resulting in the deposition of urates in the soft tissues, especially in the joint areas. Increased levels may be also found associated with leukemia, toxemia of pregnancy and severe renal impairment.

Less common are the cases of *hypouricemia* where the concentration of uric acid is below 2.0 mg/dL (0.12 mmol/L). These cases are usually secondary to cases of hepatocellular disease, renal reabsorption defect, or overtreatment with uricosuric drugs used in the treatment of hyperuricemia.

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

NOTES

- Uric acid in urine may be assayed on fresh random or timed (24-h) samples. To prevent urate precipitation specimens are brought to pH > 8 with 0.01N NaOH. Dilute urine 1:100 with distilled water before the analysis.

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