ALKALINE PHOSPHATASE BR

DGKC


Colorimetric method

KINETIC

REF KR10042 5 x 40 mL	REF KR10045 11 x 40 mL
CONTENTS R1.Reagent 5 x 32 mL R2.Reagent 5 x 8 mL	CONTENTS R1.Reagent 11 x 32 mL R2.Reagent 11 x 8 mL
For <i>in vitro</i> diagnostic use only	

SUMMARY

Alkaline phosphatase (ALP) catalyze the hydrolysis of 4 p-nitrophenylphosphate (4-NPP) with the formation of free 4-nitrophenol and inorganic phosphate, acting the alkaline buffer as a phosphate-group acceptor. The reaction is monitored kinetically at 405 nm by the rate of formation of 4-nitrophenol, proportional to the activity of ALP present in the sample.

The method follows the proposed optimised formulation of the DGKC.¹

REAGENTS

R1. ALP buffer. DEA buffer 1.25 mol/L pH 10.2, magnesium chloride 0.6 mmol/L. Biocides.

R2. ALP substrate.. 4-NPP 50 mmol/L. Biocides.

PREPARATION

The Reagents are ready-to-use.

STORAGE AND STABILITY

Store at 2-8°C. The Reagents are stable until the expiry date stated on the label. After daily use stored tightly closed and protected from light. On board these reagents are stable 3 days. Discard the reagent if the blank presents an absorbance over 1.000 at 405 nm against distilled water.

SAMPLE COLLECTION

Serum or heparinized plasma, free of hemolysis. Other anticoagulants such as EDTA, oxalate and citrate inhibit the enzyme by complexing Mg²⁺ and should not be used.

Alkaline phosphatase in serum or plasma is stable for 7 days at 2-8°C.

INTERFERENCES

- Lipemia (intralipid 20 g/L) do not interfere.
- Bilirubin (20 mg/dL) do not interfere.
- Hemoglobin (>2 g/L) may affect the results.
- Other drugs and substances may interfere^{2,3}.

INSTRUMENTATION AND MATERIALS

- KROMA analyzer.
- Laboratory equipment.
- iso-Clean Solution. Ref. CT18002.
- Multicalibrator CC/H 10x3 mL Ref. CT19750

AUTOMATED PROCEDURE

A graphic display pictures the specific sets corresponding to the technical application outlined for this test.

Any new application, to the instrument should be validated to confirm that results meet the analytical performance of the method.

It is recommended to validate periodically the instrument.

CALIBRATION

No calibration is required. Activity is calculated against theoretical factor.

CALCULATION

Patient values are calculated automatically by theoretical factor.

RESULTS

Samples with ΔA/min exceeding 0.250 at 405 nm should be diluted 1:2 with saline and assayed again. Multiply the results by 2.

If results are to be expressed as SI units apply: U/L x 0.01667 = μkat/L.

EXPECTED VALUES³

Serum, plasma

Children, up to	480 U/L (8.0 μktal/L)	25°C
Adults, up to	180 U/L (3.0 μktal/L)	25°C
Children, up to	590 U/L (9.8 μktal/L)	30°C
Adults, up to	220 U/L (3.7 μktal/L)	30°C
Children, up to	800 U/L (13.3 μktal/L)	37°C
Adults, up to	270 U/L (4.5 μktal/L)	37°C

It is recommended that each laboratory establishes its own reference range.

QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) with assayed values handled as unknowns.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

DIAGNOSTIC CHARACTERISTICS

Serum ALP measurements are of particular interest in the investigation of two groups of conditions: *bone disease* and *hepatobiliary disease*.

Among the bone diseases, the highest levels are found in Paget's disease and in patients with osteogenic bone cancer, and moderate raises in osteomalacia and rickets, the latter falling to normal on treatment with vitamin D.

Physiological bone growth elevates ALP in serum of growing children and a transient elevation may be found during healing of bone fractures.

Causes of decreased plasma ALP level are: cretinism, vitamin D deficiency and hypophosphatasia, an hereditary bone disease.

The response to the liver to any form of biliary tree obstruction is to synthesize more ALP. Intrahepatic obstruction of the bile flow by invading cancer or drugs raises serum ALP. Any drug that is hepatotoxic or induces cholestasis will greatly increase serum ALP. Well over 200 drugs have been shown to increase serum ALP in susceptible patients.⁴

PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

BIBLIOGRAPHY

1. German Society for Clinical Chemistry: Recommendations of the Enzyme Commission. Z. Klin. Chem. Klin. Biochem. 10 : 281 (1972).
2. Young, D.S. Effects of Drugs on Clinical Laboratory Tests. 4th Edition. AACC Press (1995).
3. Tietz. N.W. Clinical Guide to Laboratory Tests, 3rd Edition. W.B. Saunders Co. Philadelphia, PA. (1995).
4. Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AACC Press, 2000.

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REFERENCIAS

1. German Society for Clinical Chemistry: Recommendations of the Enzyme Commission. Z. Klin. Chem. Klin. Biochem. 10 : 281 (1972).
2. Young, D.S. Effects of Drugs on Clinical Laboratory Tests. 4th Edition. AACC Press (1995).
3. Tietz, N.W. Clinical Guide to Laboratory Tests, 3rd Edition. W.B. Saunders Co. Philadelphia, PA. (1995).
4. Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AACC Press, 2000.

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General Info

Name: Alkaline Phosph. Code: ALP Barcode: a2 Unit: U/l Decimal digits: 1

Visible: Type: Kinetic Na+ No. of rgt: 2 Multiply diluted result:

Notes:

Filters

F1: 405nm F2: not used

Incubation/Reading time [sec]

Substrate/Sample Start: Sample Start

R1,S -> R2: 36 R1,R2,S -> R3: 36 Final incub.: 120

Kinetic reading time: 192

Reference Ranges

Sample type: Serum

Patien...	Min	Max
Female	0.00	0.00
Male	0.00	0.00
Paediatric	0.00	0.00

Volumes [microlitres]

Sample: 4

R1: 160 R2: 40 R3: 0

Bottle sizes

R1: 30 ml R2: 20 ml R3: 20 ml

Kinetic/Fixed Time data

Substrate depletion: 0.5

Fit limit: 0.95

Instrument Factor (Y = aX + b)

a: 1.000 b: 0.000

Reagents

Include blank in calc.

Blank Abs (min; max): 0.3 0.9

Reagent linearity: 800

Detection limit: 9.8

Calibrations

No. of Standards: 0

Controls

C1 C2 C3

Dilutions by sample

Formula Builder

Print

Delete

Save

Name	N. of Standards	N. of rep.	Units	Stability
Alkaline Phosphatase BR	0	1	U/L	

Dil. ratio	Std value	O.D.	Reagent ...
▼			
Factor	Factor Min	Factor Max	
2770.000			

KR10042-45-1/0907