AST/GOT opt.

IFCC

UV enzymatic method


KINETIC

REF KR10072 5 x 40 mL	REF KR10075 11 x 40 mL
CONTENTS R1.Reagent 5 x 32 mL R2.Reagent 5 x 8 mL	CONTENTS R1.Reagent 11 x 32 mL R2.Reagent 11 x 8 mL
For <i>in vitro</i> diagnostic use only	

SUMMARY

Aspartate aminotransferase (AST/GOT) catalyzes the transfer of the amino group from aspartate to oxoglutarate with the formation of glutamate and oxalacetate. The latter is reduced to malate by malate dehydrogenase (MDH) in the presence of reduced nicotinamide adenine dinucleotide (NADH).

The reaction is monitored kinetically at 340 nm by the rate of decrease in absorbance resulting from the oxidation of NADH to NAD⁺, proportional to the activity of AST present in the sample.

This test has been formulated according the standardized method described by IFCC. Clin Chem Lab Med 2002; 40(7) : 718-724².

REAGENTS

R1. AST substrate. TRIS buffer 121 mmol/L pH 7.8, L-aspartate 362 mmol/L, malate dehydrogenase > 460 U/L, lactate dehydrogenase > 600 U/L.

R2. AST coenzyme. NADH 1.3 mmol/L, 2-oxoglutarate 75 mmol/L. Biocides.

PREPARATION

The Reagents are ready-to-use.

STORAGE AND STABILITY

Store at 2-8°C. The Reagents are stable until the expiry date stated on the label. After daily use stored tightly closed and protected from light.

On board the reagents are stable 30 days.

Discard the reagent if the blank presents an absorbance lower 1.300 at 340 nm against distilled water.

SAMPLE COLLECTION

Serum and EDTA or heparinized plasma free of hemolysis.

AST is stable in serum or plasma 24 hours at room temperature and for 1 week at 2-8°C.

INTERFERENCES

- Lipemia (intralipid >15 g/L) does not interfere.
- Bilirubin (>30 mg/dL) does not interfere.
- Hemoglobin (>10 g/dL) does not interfere.
- Other drugs and substances may interfere².

INSTRUMENTATION AND MATERIALS

- KROMA analyzer.
- Laboratory equipment.
- iso-Clean Solution. Ref. CT18002.
- Multicalibrator CC/H 10x3 mL Ref. CT19750

AUTOMATED PROCEDURE

A graphic display pictures the specific sets corresponding to the technical application outlined for this test.

Any new application, to the instrument should be validated to confirm that results meet the analytical performance of the method.

It is recommended to validate periodically the instrument.

CALIBRATION

No calibration is required. Activity is calculated against theoretical factor.

CALCULATION

Patient values are calculated automatically by theoretical factor.

RESULTS

Samples with ΔA/min exceeding 0.160 at 340 nm should be diluted 1:10 with saline and assayed again. Multiply the results by 10.

If results are to be expressed as SI units apply: U/L x 0.01667 = μkat/L

EXPECTED VALUES³

Serum, plasma

Adults	37°C	up to 40 U/L (0.67 μkat/L)
	30°C	up to 25 U/L (0.42 μkat/L)

Levels approximately twice the adult level are seen in neonates and infants; these decline to adult level by approximately 6 months of age.

It is recommended that each laboratory establishes its own reference range.

QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) with assayed values handled as unknowns.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

DIAGNOSTIC CHARACTERISTICS

The group of enzymes called transaminase exist in tissues of many organs. Necrotic activity in these organs causes a release of abnormal quantities of enzyme into the blood where they are measured.

Since heart tissue is rich in AST increased serum levels appear in patients after myocardial infarction, as well as in patients with muscle disease, muscular dystrophy and dermatomyositis. The liver is specially rich in ALT, being this enzyme measurement used primarily as a test for infectious and toxic hepatitis, although high levels of both ALT and AST may also be found in cases of liver cell damage and acute pancreatitis, suggesting that the obstruction of the biliary tree by the edematous pancreas and the presence of associate hepatic disease may contribute to elevated AST levels in these patients. Slight or moderate elevations of AST and ALT activities may be observed after intake of alcohol and after administration of various drugs, such as salicylates, opiates and ampicillin. Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

BIBLIOGRAPHY

1. Bergmeyer, H.V., Hørdler, M., Rej, R. Approved recommendation (1985) on IFCC methods for the measurement of catalytic concentration of enzymes. Part 2. IFCC method for aspartate aminotransferase, J. Clin. Chem. Clin. Biochem. 24 : 497 (1986).
2. Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AACC Press, 2000.
3. Tietz. N.W. Clinical Guide to Laboratory Tests, 3rd Edition. W.B. Saunders Co. Philadelphia, PA. (1995).

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REFERENCIAS

1. Bergmeyer, H.V., Hørdler, M., Rej, R. Approved Recommendation (1985) on IFCC Methods for the Measurement of Catalytic Concentration of Enzymes. Part 2. IFCC Method for Aspartate Aminotransferase, J. Clin. Chem. Clin. Biochem. 24 : 497 (1986).
2. Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AACC Press, 2000.
3. Tietz. N.W. Clinical Guide to Laboratory Tests, 3rd Edition. W.B. Saunders Co. Philadelphia, PA. (1995).

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General Info

Name: AST / GOT Code: GOT Barcode: a5 Unit: U/l Decimal digits: 2

Visible: Type: Kinetic Na+ No. of rgt: 2 Multiply diluted result:

Notes:

Filters

F1: 340nm F2: not used

Volumes [microlitres]

Sample: 11

R1: 160 R2: 40 R3: 0

Bottle sizes

R1: 30 ml R2: 20 ml R3: 20 ml

Incubation/Reading time [sec]

Substrate/Sample Start: Sample Start

R1,S -> R2: 36 R1,R2,S -> R3: 36 Final incub.: 60

Kinetic reading time: 192

Kinetic/Fixed Time data

Substrate depletion: 0.5 Fit limit: 0.95

Reference Ranges

Sample type: Serum

Patient...	Min	Max
Female	0.00	0.00
Male	0.00	0.00
Paediatric	0.00	0.00

Instrument Factor (Y = aX + b)

a: 1.000 b: 0.000

Controls

C1 C2 C3

Reagents

Include blank in calc.

Blank Abs (min; max): 0.8 2

Reagent linearity: 500

Detection limit: 2.4

Calibrations

No. of Standards: 0

Buttons: Formula Builder, Print, Delete, Save, Dilutions by sample

Name	N. of Standards	N. of rep.	Units	Stability
AST/GOT BR	0	1	U/L	

Dil.ratio	Std value	O.D.	Reagent ...

Factor	Factor Min	Factor Max
-3333.000		

KR10072-75-1/0907