## BILIRUBIN

TOTAL AND DIRECT  
*Colorimetric method*

END POINT

REF KR10082 RD. 5 x 40 mL CONTENTS RD.Reagent 5 x 32 mL RN.Reagent 5 x 8 mL	REF KR10085 RD. 11 x 40 mL CONTENTS RD.Reagent 11 x 32 mL RN.Reagent 11 x 8 mL	REF KR10092 RT. 5 x 40 mL CONTENTS RT.Reagent 5 x 32 mL RN.Reagent 5 x 8 mL	REF KR10095 RT. 11 x 40 mL CONTENTS RT.Reagent 11 x 32 mL RN.Reagent 11 x 8 mL
For <i>in vitro</i> diagnostic use only			

### SUMMARY

Bilirubin is converted to coloured azobilirubin by diazotized sulfanilic acid and is measured photometrically. Of the two bilirubin fractions in serum –bilirubin-glucuronide and free bilirubin which is bound to albumin– only the former reacts directly, while free albumin reacts after being displaced from protein by an accelerator. The difference of two measurements total bilirubin (with accelerator) and direct bilirubin (without accelerator) enables to calculate indirect bilirubin. The terms «direct» and «indirect» bilirubin refers exclusively to the reaction characteristics in the presence or absence of an accelerator or solubilizer and are only approximate equivalents of the two bilirubin fractions.<sup>1,2</sup>

### REAGENTS

**RT.** Sulfanilic acid 29 mmol/L, HCl 0.24 mol/L, Duposol® 3% (w/v).

**RD.** Sulfanilic acid 29 mmol/L, HCl 0.24 mol/L.

**RN.** Sodium nitrite 11.6 mmol/L.

### PREPARATION

The Reagents are ready-to-use.

### STORAGE AND STABILITY

 Store at 2-8°C. The Reagents are stable until the expiry date stated on the label. On board the reagents are stable 30 days. After daily use stored tightly closed and protected from light.

#### Discard if appear signs of deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 540 nm > 0.050 in 1cm cuvette.
- Reagent **RN** if it develops a yellow coloration.

### SAMPLE COLLECTION

Fresh hemolysis-free serum, EDTA or heparinized plasma. Store in the dark until use.

Samples can be frozen at -15°C or below in which case bilirubin is stable for 2 months.

### INTERFERENCES

- Lipemia (intralipid < 5 g/L) does not interfere.
- Direct Bilirubin (Hemoglobin 2 g/L) may affect the results.
- Total Bilirubin (Hemoglobin 16 g/L) does not interfere.
- Other drugs and substances may interfere<sup>3</sup>.
- Lipemic samples interfere with the assay. The interference can be corrected by preparing a sample blank before applying the general formula of calculation.

### INSTRUMENTATION AND MATERIALS

- KROMA analyzer.
- Laboratory equipment.
- iso-Clean Solution. Ref. CT18002.
- Multicalibrator CC/H 10x3 mL Ref. CT19750

### AUTOMATED PROCEDURE

A graphic display pictures the specific sets corresponding to the technical application outlined for this test.

Any new application, to the instrument should be validated to confirm that results meet the analytical performance of the method.

It is recommended to validate periodically the instrument.

### CALIBRATION

Recalibrate weekly (Bilirubin T) or every 2 days (Bilirubin D), when a new lot of reagent is used, when control recovery falls out of the expected range or when adjustments are made to the instrument. Two point calibration is recommended (S1: NaCl 9 g/L and S2: Calibrator). A reagent blank should be run daily before sample analysis.

### RESULTS

Samples with concentrations higher than 20 mg/dL should be diluted 1:2 with saline and assayed again. Multiply results by 2.

If results are to be expressed as SI units apply: mg/dL x 17.1 = µmol/L

### EXPECTED VALUES <sup>4</sup>

ADULTS Total Up to 1.0 mg/dL	Direct Up to 0.2 mg/dL
------------------------------	------------------------

### NEWBORNS (TOTAL BILIRUBIN)

Age	Premature	Full- term
Up to 24 h	1.0-6.0 mg/dL	2.0-6.0 mg/dL
Up to 48 h	6.0-8.0 mg/dL	6.0-7.0 mg/dL
3-5 days	10.0-15.0 mg/dL	4.0-12.0 mg/dL

It is recommended that each laboratory establishes its own reference range.

### QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) with assayed values handled as unknowns.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

### DIAGNOSTIC CHARACTERISTICS

*Hyperbilirubinemia* (an abnormal elevation of bilirubin, whether conjugated or unconjugated) in plasma is an indication of a disturbance in bilirubin metabolism. This condition is caused either by an overproduction of bilirubin or by an impairment in the metabolic pathway.

The increase in bilirubin production is usually caused by a rapid destruction of erythrocytes, resulting from blood diseases such as haemolytic anemia.

In newborns the increase in bilirubin may be caused by Rh, ABO, or other blood group incompatibility, by sepsis, hepatic immaturity, or by a variety of hereditary defects in bilirubin conjugation.

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

### PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

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*In vitro diagnostic reagent for the determination of total bilirubin in serum or plasma.*

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*The reagents **RT** and **RN** are ready to use.*

**CALIBRADOR**  
*STANDARD*

**Linear Multicalibrator CC/H, Code CT19750**

**CONTROL CALIDAD**  
*QUALITY CONTROL*

**Linear Human Multisera Normal CT19800**

**Linear Human Multisera Abnormal CT19850**

**DETERMINACIONES**  
*NUMBER OF TESTS*

**1000 tests/kit, Code KR10082**

**2.200 tests/kit, Code KR10085**

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*(dead volume is not taken in consideration).*

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*Without automatic post-dilution: Up to 10 mg/dL*

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*For additional information read reagent packaging insert.*

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*The reported application is given only as a guideline. Any application to an instrument should be validated to demonstrate that results meet the performance characteristics of the method.*

**General Info**

Name:  Code:  Barcode:  Unit:  Decimal digits:

Visible Type:   No. of rgt:  Multiply diluted result:

Notes:

**Filters**

F1:  F2:

**Incubation/Reading time [sec]**

Substrate/Sample Start

Sample Start:

R1,S -> R2    R1,R2,S -> R3    Final incub.

**Reference Ranges**

Sample type:

Patien...	Min	Max
Female	0.30	1.00
Male	0.30	0.80
Paediatric	0.60	1.90

**Volumes [microlitres]**

Sample:

R1:  R1:

R2:  R2:

R3:  R3:

**Kinetic/Fixed Time data**

**Instrument Factor (Y = aX + b)**

a:  b:

**Reagents**

Include blank in calc.

Blank Abs (min; max):

Reagent linearity:

Detection limit:

**Calibrations**

No. of Standards:

Piecewise Linear  Calib. Curve

**Controls**

C1     C2     C3

Dilutions by sample

Formula Builder

Print

Delete

Save

Name	N. of Standards	N. of rep.	Units	Stability
Direct Bilirubin	2	1	mg/dL	

Dil.ratio	Std value	O.D.	Reagent ...
<input type="text"/>	*	<input type="text"/>	<input type="text"/>
Factor	Factor Min	Factor Max	
<input type="text"/>	<input type="text"/>	<input type="text"/>	

\* entered by the user

KR10082-85-1/0907