CHOLESTEROL MR

TOTAL


Enzymatic colorimetric method

ENDPOINT

REF KR10142 8 x 30 mL	REF KR10145 18 x 30 mL
CONTENTS R1.Reagent 8 x 30 mL	CONTENTS R1.Reagent 18 x 30 mL
For <i>in vitro</i> diagnostic use only	

SUMMARY

This method for the measurement of total cholesterol^{1,2} in serum involves the use of three enzymes: cholesterol esterase (CE), cholesterol oxidase (CO) and peroxidase (POD). In the presence of the former the mixture of phenol and 4-aminoantipyrine (4-AA) are condensed by hydrogen peroxide to form a quinoneimine dye proportional to the concentration of cholesterol in the sample.

REAGENT

R1. PIPES 200 mmol/L pH 7.0, sodium cholate 1 mmol/L, cholesterol esterase > 250 U/L, cholesterol oxidase > 250 U/L, peroxidase > 1 KU/L, 4-aminoantipyrine 0.33 mmol/L, phenol 4 mmol/L, non-ionic tensioactives 2 g/L (w/v). Biocides.

PREPARATION

The Monoreagent is ready-to-use.

STORAGE AND STABILITY

Store at 2-8°C. The Reagent is stable until the expiry date stated on the label. On board the reagent is stable 30 days. After daily use stored tightly closed and protected from light. Discard the reagent if presents an absorbance over 0.200 at 500 nm against distilled water.

SAMPLE COLLECTION

Serum, EDTA or heparinized plasma free of hemolysis. Cholesterol in serum or plasma is stable up to 5 days at 2-8°C and for a few months at -20°C.

INTERFERENCES

- Lipemia (intralipid 5 g/L) may affect the results.
- Bilirubin (40 mg/dL) does not interfere.
- Hemoglobin (> 1 g/L) may affect the results.
- Other drugs and substances may interfere³.

INSTRUMENTATION AND MATERIALS

- KROMA analyzer.
- Laboratory equipment.
- iso-Clean Solution. Ref. CT18002.
- Multicalibrator CC/H 10x3 mL Ref. CT19750

AUTOMATED PROCEDURE

A graphic display pictures the specific sets corresponding to the technical application outlined for this test. Any new application, to the instrument should be validated to confirm that results meet the analytical performance of the method. It is recommended to validate periodically the instrument.

CALIBRATION

Recalibrate weekly, when a new lot of reagent is used, when control recovery falls out of the expected range or when adjustments are made to the instrument. A reagent blank should be run daily before sample analysis.

RESULTS

Samples with concentrations higher than 600 mg/dL should be diluted 1:2 with saline and assayed again. Multiply the results by 2. If results are to be expressed as SI units apply: mg/dL x 0.0259 = mmol/L

EXPECTED VALUES⁴

Updated clinical values of total cholesterol used to classify risk groups.

Total Cholesterol	Risk Classification
< 200 mg/dL (< 5.18 mmol/L)	Desirable
200-239 mg/dL (5.18-6.2 mmol/L)	Borderline high
> 240 mg/dL (> 6.2 mmol/L)	High

It is recommended that each laboratory establishes its own reference range.

QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) with assayed values handled as unknowns. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

DIAGNOSTIC CHARACTERISTICS

The knowledge of the plasma level of lipids (cholesterol and triglycerides) together with lipoproteins of high and low density (HDL and LDL) aids in the detection of many conditions bound to metabolic disorders of high risk. The imbalance in the level of lipoproteins in plasma leads to *hyperlipoproteinemias*, a group of disorders that affects lipid levels in serum, causing coronary heart disease (CHD) and atherosclerosis, conditions in which the cholesterol levels are important tools in their diagnosis and classification.

Jaundice of the obstructive type usually is accompanied by an elevated total serum cholesterol with a normal ester fraction. Diabetes, hypothyroidism, and certain types of kidney disease are other disorders that may exhibit the same cholesterol disturbance.

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

BIBLIOGRAPHY

1. Allain, C.C., Poon, L.S., Clau, C.S.G, Richmond, W and Fu, P.D. Clin. Chem. 20 : 470 (1974).
2. Tamaoku, K., Murao, Y. and Akiura, K. Analytica Chimica. Acta. 136 : 121 (1982).
3. Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AACC Press, 2000.
4. SPECIAL REPORT. Executive Summary of the Third Report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). JAMA. 285 : 2486 (2001).

Further hints

- SPECIAL REPORT (ATP III) available at:
<http://www.nhlbi.nih.gov>
- An autoevaluation about the risk of heart disease is available at :
time.com/cholesterol

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REFERENCIAS

1. Allain, C.C., Poon, L.S., Clau, C.S.G, Richmond, W y Fu, P.D. Clin. Chem. 20 : 470 (1974).
2. Tamaoku, K., Murao, Y. y Akiura, K. Analytica Chimica. Acta. 136 : 121 (1982).
3. Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AACC Press, 2000.
4. SPECIAL REPORT. Executive Summary of the Third Report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). JAMA. 285 : 2486 (2001).

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General Info

Name: Cholesterol Code: CHOL Barcode: a9 Unit: mg/dl Decimal digits: 2

Visible: Type: Endpoint Na+ No. of rgt: 1 Multiply diluted result:

Notes:

Filters

F1: 505nm F2: not used

Volumes [microlitres]

Sample: 2

R1: 200 R2: 0 R3: 0

Bottle sizes

R1: 30 ml R2: 20 ml R3: 20 ml

Incubation/Reading time [sec]

Substrate/Sample Start: Substrate Start

R1,S -> R2: 36 R1,R2,S -> R3: 36 Final incub.: 300

Reference Ranges

Sample type: Serum

Patien...	Min	Max
Female	0.00	0.00
Male	0.00	0.00
Paediatric	0.00	0.00

Instrument Factor (Y = aX + b)

a: 1.000 b: 0.000

Controls

C1 C2 C3

Reagents

Include blank in calc.

Blank Abs (min; max): 0.001 0.3

Reagent linearity: 600

Detection limit: 2.8

Calibrations

No. of Standards: 2

Piecewise Linear Calib. Curve

Buttons: Formula Builder, Print, Delete, Save

Name	N. of Standards	N. of rep.	Units	Stability
Cholesterol MR	2	1	mg/dL	

Dil.ratio	Std value	O.D.	Reagent ...
	*		

Factor	Factor Min	Factor Max

* entered by the user

KR10142-45-1/0907