# CREATINE KINASE MB

CK-MB INHIBITED

*Immunological UV method*


KINETIC

<b>REF KR10170</b> 1 x 40 mL	<b>REF KR10172</b> 5 x 40 mL
<b>CONTENTS</b> R1.Reagent 1 x 32 mL R2.Reagent 1 x 8 mL	<b>CONTENTS</b> R1.Reagent 5 x 32 mL R2.Reagent 5 x 8 mL
For <i>in vitro</i> diagnostic use only	

## SUMMARY

The major CK activity in normal serum is due to CK-MM and CK-MB isoenzymes, found mainly in skeletal and heart muscle. CK-BB is usually present in serum at low concentration.

Both, creatine kinase (CK) enzymes are dimers formed by the association of two sub-units from muscle (M) and nerve cells (B). The immunoinhibition from an specific antibody of both, MM sub-units and the single M sub-unit of CK-MB, allows the determination of the B sub-unit. The CK-B activity corresponding to half of CK-MB is measured by the increasing rate of absorbance resulting from the following coupled reactions.<sup>1,2</sup>

## REAGENTS

**R1. Buffer/Glucose/NAC.** Imidazol buffer 100 mmol/L pH 6.7, glucose 20 mmol/L, NAC 20 mmol/L, magnesium acetate 10 mmol/L, NADP 2.5 mmol/L, HK  $\geq 4$  KU/L, EDTA 2 mmol/L.

**R2. Substrate/Coenzymes.** CP 30 mmol/L, AMP 5 mmol/L, ADP 2 mmol/L, di(adenosine-5') pentaphosphate 10  $\mu$ mol/L, G6P-DH  $\geq 1.5$  KU/L. Sufficient CK-M human antibody to inhibit  $\geq 3000$  U/L of CK-MM at 37°C.

## PREPARATION

**Working reagent.** Mix 4 mL of R1 + 1 mL of R2.

## STORAGE AND STABILITY

Store at 2-8°C. The Reagents are stable until the expiry date stated on the label. After daily use stored tightly closed and protected from light. On board are stable 15 days. Discard the reagents if presents an absorbance over 0.400 at 340 nm against distilled water.

## SAMPLE COLLECTION

Serum. Stable 8 days at 2-8°C or 1 month at -20°C. Chill the samples as rapidly as possible after collection.

Moderately or severely hemolyzed specimens are unsatisfactory for testing, as well as plasmas containing EDTA, heparin, citrate, or fluoride since they may produce unpredictable reaction rates.

## INTERFERENCES

- Lipids (intralipid >1 g/L) may affect the results.
- Bilirubin (< 20 mg/dL) does not interfere.
- Hemoglobin (>8 g/L) may affect the results.
- The presence of above normal concentrations of CK-BB or adenilatekinasa, and of macro or mitochondrial CK may affect the results<sup>6</sup>.
- Other drugs and substances may interfere<sup>8-9</sup>.

## INSTRUMENTATION AND MATERIALS

- KROMA analyzer.
- Laboratory equipment.
- iso-Clean Solution. Ref. CT18002.
- Multicalibrator CC/H 10x3 mL Ref. CT19750

## AUTOMATED PROCEDURE

A graphic display pictures the specific sets corresponding to the technical application outlined for this test. Any new application, to the instrument

## CALIBRATION

No calibration is required. Activity is calculated against theoretical factor.

## CALCULATION

Patient values are calculated automatically by theoretical factor.

## RESULTS

*Percentage of CK-MB activity:*

$$\frac{\text{CK-MB Activity}}{\text{Total CK Activity}} \times 100 = \% \text{ CK-MB activity}$$

## EXPECTED VALUES

The normal range of CK-MB activity in the serum of adults is considered to be 2.0 U/L to 19.5 U/L at 37°C. Newborns, infants, and children have higher serum CK-MB values than adults.

A ratio between CK-MB and total CK activities above 4% should be considered suspicious, and above 10% consistent with acute myocardial infarction.<sup>4</sup> When the three conditions stated below are met the probability of myocardial injury is very high.<sup>5</sup>

Condition	Test temperature (37°)	
	1	CK Men
	CK Women	> 167 U/L (2.78 $\mu$ kat/L)
2	CK-MB	> 24 U/L (0.40 $\mu$ kat/L)
3	% CK-MB	6-25%

It is recommended that each laboratory establishes its own reference range.

## QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) with assayed values handled as unknowns. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

## DIAGNOSTIC CHARACTERISTICS

The major medical application of CK-MB assays is in adults for the diagnosis of acute myocardial infarction (AMI) or, the differentiation of myocardial injury from skeletal muscle damage. Creatine kinase-MB is considered one of the best laboratory indicators of AMI. CK-MB determination within the proper time frame after infarction is most critical being the useful interval ("window") from about 10 to 24 hour after the infarction. Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

## PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

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*In vitro diagnostic reagent for the determination of creatine kinase MB (CK-MB) in serum.*

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*Factor included in the application a 37°C.*

**CONTROL CALIDAD**  
*QUALITY CONTROL*

**Linear Human Multisera Normal CT19800**

**Linear Human Multisera Abnormal CT19850**

**DETERMINACIONES**  
*NUMBER OF TESTS*

**200 tests/kit, Code KR10170**

**1000 tests/kit, Code KR10172**

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*(dead volume is not taken in consideration).*

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*Without automatic post-dilution: Up to 330 U/L*

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*For additional information read reagent packaging insert.*

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*The reported application is given only as a guideline. Any application to an instrument should be validated to demonstrate that results meet the performance characteristics of the method.*

**General Info**

Name  Code  Barcode  Unit  Decimal digits

Visible Type   No. of rgt  Multiply diluted result

Notes

**Filters**

F1  F2

**Incubation/Reading time [sec]**

Substrate/Sample Start

R1,S -> R2    R1,R2,S -> R3    Final incub.

Kinetic reading time

**Reference Ranges**

Sample type

Patien...	Min	Max
Female	0.00	0.00
Male	0.00	0.00
Paediatric	0.00	0.00

**Volumes [microlitres]**

Sample

R1  R2  R3

**Bottle sizes**

R1  R2  R3

**Kinetic/Fixed Time data**

Substrate depletion

Fit limit

**Instrument Factor (Y = aX + b)**

a  b

**Reagents**

Include blank in calc.

Blank Abs (min; max)

Reagent linearity

Detection limit

**Calibrations**

No. of Standards

**Controls**

C1     C2     C3

Dilutions by sample

Name	N. of Standards	N. of rep.	Units	Stability
Creatine Kinase MB	0	1	U/L	

Dil.ratio	Std value	O.D.	Reagent ...
▼			
Factor	Factor Min	Factor Max	
8360.000	0.000	0.000	

KR10170-72-1/0907