HDL-CHOLESTEROL

DIRECT

Enzymatic colorimetric method

FIXED TIME

REF KR10242 5 x 40 mL CONTENTS R1.Reagent 5 x 30 mL R2.Reagent 5 x 10 mL	REF KR10245 11 x 40 mL CONTENTS R1.Reagent 11 x 30 mL R2.Reagent 11 x 10 mL
For <i>in vitro</i> diagnostic use only	

SUMMARY

This direct method^{1,2} for quantifying cholesterol in high-density lipoproteins (HDL) is a homogeneous enzymatic test in which the differential precipitation and further sedimentation of the rest of lipoproteins and quilomicrons is avoided.

The procedure comprises two steps. In the first step cholesterol in lipoproteins other than HDL in the test sample are decomposed by the simultaneous action of cholesterol esterase (CE) and cholesterol oxidase (CO) at pH 7.0, giving as end products cholestenone and hydrogen peroxide, the latter being decomposed to water and oxygen by catalase.

In the second step a surfactant which specifically acts on HDL is added to the reaction product of the first step being the remaining cholesterol quantified by a Trinder's type reaction in which the aniline derivate, HDAOS*, and 4-aminoantipyrine (4-AA) as a coupling reagent are condensed by the H₂O₂ in presence of peroxidase (POD) to form a red quinoneimine dye proportional to the concentration of HDL-cholesterol present in the sample.

* *N*-(2-hydroxy-3-sulfopropyl)-3,5-dimethoxyaniline

REAGENTS

R1. Enzyme. GOOD's buffer 100 mmol/L pH 7.0, MgCl₂ 18 mmol/L, CE 800 U/L, CO 500 U/L, catalase 100 KU/L, HDAOS 0.7 mmol/L.

R2. POD/4-AA reagent. POD 4 KU/L, 4-AA 4 mmol/L, N₃Na < 0.1%, specific surfactants < 1.5 % (v/v).

PREPARATION

The Reagents are ready-to-use.

STORAGE AND STABILITY

Store at 2-8°C. The Reagents are stable until the expiry date stated on the label. After daily use stored tightly closed and protected from light. On board the reagents are stable 30 days. Discard if appear signs of deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 600 nm < 0.080 in 1cm cuvette.

SAMPLE COLLECTION

Serum, EDTA or heparinized plasma obtained by the patient after an overnight fast. Remove from cells within 3 hours of venipuncture. Samples may be kept at 4-8°C for 2 weeks or at -20°C for 3 months.

INTERFERENCES

- Lipemia (Triglycerides < 12 g/L) does not interfere.
- Bilirubin (< 30 mg/dL) does not interfere.
- Hemoglobin (< 5 mg/dL) does not interfere.
- Ascorbate (< 50 mg/dL) does not interfere.
- Other drugs and substances may interfere⁶⁻⁷.

INSTRUMENTATION AND MATERIALS

- KROMA analyzer.
- Laboratory equipment.
- iso-Clean Solution. Ref. CT18002.
- Multicalibrator CC/H 10x3 mL Ref. CT19750

AUTOMATED PROCEDURE

A graphic display pictures the specific sets corresponding to the technical application outlined for this test. Any new application, to the instrument should be validated to confirm that results meet the analytical performance of the method. It is recommended to validate periodically the instrument.

CALIBRATION

Recalibrate weekly, when a new lot of reagent is used, when control recovery falls out of the expected range or when adjustments are made to the instrument. Two point calibration is recommended (S1: NaCl 9 g/L and S2: Calibrator). A reagent blank should be run daily before sample analysis.

RESULTS

If results are to be expressed as SI units apply:
mg/dL x 0.0259 = mmol/L

EXPECTED VALUES³⁻⁵

Clinical values of HDL-Cholesterol used to classify risk groups.

Cholesterol from lipoproteins of high density		RISK
Men	> 55 mg/dL (> 1.42 mmol/L)	Low
	40-55 mg/dL (0.90-1.42 mmol/L)	Moderate
	< 40 mg/dL (< 1.04 mmol/L)	High
Women	> 65 mg/dL (> 1.68 mmol/L)	Low
	45-65 mg/dL (1.16-1.68 mmol/L)	Moderate
	< 45 mg/dL (< 1.16 mmol/L)	High

QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) with assayed values handled as unknowns. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

DIAGNOSTIC CHARACTERISTICS

Low HDL-cholesterol is a strong independent predictor of coronary heart disease. In ATP III⁴, low HDL cholesterol is defined categorically as a level < 40 mg/dL (1.04 mmol/L), a change from the level of < 35 mg/dL in ATPII (1993).

Low HDL cholesterol is used as a risk factor to estimate 10-year risk for coronary heart disease, having several causes: elevated triglycerides, overweight and obesity, physical inactivity, and type 2 diabetes.

Other causes are, cigarette smoking, very high carbohydrate intakes (> 60% of calories), and certain drugs as anabolic steroids and progestational agents.

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

BIBLIOGRAPHY

1. US Patent No.: 6,479,249 B2 (Nov 12,2002).
2. Sachiko Izawa. J.Med. and Pharm. Sci. 37 : 1325 (1997).
3. Warnick, G. y Wood, P.D., Clin. Chem. 41 : 1427 (1995).
4. SPECIAL REPORT. Executive Summary of the Third Report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). JAMA. 285 : 2486 (2001).
5. Tietz. N.W. Fundamentals of Clinical Chemistry,p.940. W.B. Saunders Co. Philadelphia, PA. (1987).
6. Friedman and Young. Effects of disease on clinical laboratory tests. 5th ed. AACC (Press 2000).
7. Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AACC Press, 2000.

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REFERENCIAS

1. US Patent No.: 6,479,249 B2 (Nov 12,2002).
2. Sachiko Izawa. J.Med. y Pharm. Sci. 37 : 1325 (1997).
3. Warnick, G. y Wood, P.D., Clin. Chem. 41 : 1427 (1995).
4. SPECIAL REPORT. Executive Summary of the Third Report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). JAMA. 285 : 2486 (2001).
5. Tietz. N.W. Fundamentals of Clinical Chemistry, p.940. W.B. Saunders Co. Philadelphia, PA. (1987).
6. Friedman y Young. Effects of disease on clinical laboratory tests. 5th ed. AACC (Press 2000).
7. Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AACC Press, 2000.

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General Info

Name: HDL Chol. Dir. Code: HDL Barcode: b6 Unit: mg/dl Decimal digits: 1

Visible: Type: Endpoint Na+ No. of rgt: 2 Multiply diluted result:

Notes:

Filters

F1: 578nm F2: not used

Volumes [microlitres]

Sample: 3

R1: 150 R2: 50 R3: 0

Bottle sizes

R1: 30 ml R2: 20 ml R3: 20 ml

Reagents

Include blank in calc.

Blank Abs (min; max): 0.001 1

Reagent linearity: 150

Detection limit: 3

Incubation/Reading time [sec]

Substrate/Sample Start: Substrate Start

R1,S -> R2: 180 R1,R2,S -> R3: 36 Final incub.: 420

Kinetic/Fixed Time data

Calibrations

No. of Standards: 2

Piecewise Linear Calib. Curve

Reference Ranges

Sample type: Serum

Patien...	Min	Max
Female	0.00	0.00
Male	0.00	0.00
Paediatric	0.00	0.00

Instrument Factor (Y = aX + b)

a: 1.000 b: 0.000

Controls

C1: C2: C3:

Dilutions by sample

Formula Builder Print Delete Save

Name	N. of Standards	N. of rep.	Units	Stability
HDL-Cholesterol Direct	2	1	mg/dL	

Dil.ratio	Std value	O.D.	Reagent ...
	*		

Factor	Factor Min	Factor Max

* entered by the user

KR10242-45-1/0907