# LDH BR

SFBC

UV enzymatic method


KINETIC

<b>REF KR10282</b> 5 x 40 mL <b>CONTENTS</b> R1. Reagent 5 x 32 mL R2. Reagent 5 x 8 mL
For <i>in vitro</i> diagnostic use only

## SUMMARY

Lactate dehydrogenase (LD/LDH) catalyzes the reduction of pyruvate to lactate (P-L) in the presence of reduced nicotinamide adenine dinucleotide (NADH) at pH 7.5.

The reaction is monitored kinetically at 340 nm by the rate of decrease in absorbance resulting from the oxidation of NADH to NAD<sup>+</sup> proportional to the activity of LD present in the sample.

The method follows the proposed optimised formulation of the SFBC.<sup>1</sup>

## REAGENTS

**R1. LDH substrate.** TRIS buffer 100 mmol/L pH 7.5, pyruvate 2.75 mmol/L, sodium chloride 222 mmol/L.

**R2. LDH coenzyme.** NADH 1.55 mmol/L.

## PREPARATION

The Reagents are ready-to-use.

## STORAGE AND STABILITY

Store at 2-8°C. The Reagents are stable until the expiry date stated on the label. After daily use stored tightly closed and protected from light. On board the reagents are stable 30 days.

Discard the reagent if the blank presents an absorbance lower 1.000 at 340 nm against distilled water.

## SAMPLE COLLECTION

Serum free of hemolysis separated from the cells as soon as possible after collection. The use of heparin and citrate as anticoagulants have been reported to falsely elevate LD activity.

Freezing results in loss of activity of the enzyme.

## INTERFERENCES

- Lipemia (intralipid > 20 g/L) does not interfere.
- Bilirubin (> 40 mg/dL) does not interfere.
- Hemoglobin (> 4 g/dL) may affect the results.
- Other drugs and substances may interfere<sup>3</sup>.

## INSTRUMENTATION AND MATERIALS

- KROMA analyzer.
- Laboratory equipment.
- Cleaning solutions Ref. KR18100, KR18200, KR18300, KR18400.

## AUTOMATED PROCEDURE

A graphic display pictures the specific sets corresponding to the technical application outlined for this test.

Any new application, to the instrument should be validated to confirm that results meet the analytical performance of the method.

It is recommended to validate periodically the instrument.

## CALIBRATION

No calibration is required. Activity is calculated against theoretical factor.

## CALCULATION

Patient values are calculated automatically by theoretical factor.

## RESULTS

Samples with ΔA/min exceeding 0.150 at 340 nm should be diluted 1:10 with saline and assayed again. Multiply the results by 10.

If results are to be expressed as SI units apply: U/L x 0.01667 = μkat/L

## EXPECTED VALUES<sup>1</sup>

Serum

Temperature	37°C	30°C
Adults	207-414 U/L (3.40-6.80 μKat/L)	140-280 U/L (2.30-4.70 μKat/L)

It is recommended that each laboratory establishes its own reference range.

## QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) Ref. CT19800-CT19850 with assayed values handled as unknowns.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

## DIAGNOSTIC CHARACTERISTICS<sup>2,5</sup>

The enzyme activity found in circulation is a mixture of five isoenzymes. Each organ has a characteristic isoenzyme profile. Leakage of these isoenzymes from a diseased organ results in the elevation of total serum LD.

*Increased* levels become evident 8-12 hours after myocardial infarction reaching a maximum 4-5 days later. Elevated values in serum are also encountered in cases of pulmonary embolism and in about one third of patients with renal disease, specially those with pyelonephritis or tubular necrosis. In toxic hepatitis with jaundice, Hodgkin's disease and abdominal and lung cancers elevations are specially high.

*Moderate* increases are also observed in cases of liver disease, megaloblastic and pernicious anemia and progressive muscular dystrophy.

*Decreases* are not important clinically.

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

## PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

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*In vitro diagnostic reagent for the determination of lactate dehydrogenase (LDH, P→L) in serum.*

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*Factor included in the application a 37°C.*

### CONTROL CALIDAD *QUALITY CONTROL*

**Linear Human Multisera Normal CT19800**

**Linear Human Multisera Abnormal CT19850**

### DETERMINACIONES *NUMBER OF TESTS*

**1000 tests/kit, Code KR10282**

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*(dead volume is not taken in consideration).*

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*Without automatic post-dilution: Up to 1500 U/L*

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*For additional information read reagent packaging insert.*

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*The reported application is given only as a guideline. Any application to an instrument should be validated to demonstrate that results meet the performance characteristics of the method.*

**General Info**

Name: LDH Code: LDH Barcode: b9 Unit: U/l Decimal digits: 1

Visible:  Type: Kinetic Na+ No. of rgt: 2 Multiply diluted result:

Notes:

---

**Filters**

F1: 340nm F2: not used

**Volumes [microlitres]**

Sample: 5

R1: 160 R2: 40 R3: 0

**Bottle sizes**

R1: 30 ml R2: 20 ml R3: 20 ml

**Reagents**

Include blank in calc.

Blank Abs (min; max): 0.8 2

Reagent linearity: 1500

Detection limit: 19

**Incubation/Reading time [sec]**

Substrate/Sample Start: Substrate Start

R1,S -> R2: 36 R1,R2,S -> R3: 300 Final incub.: 36

Kinetic reading time: 192

**Kinetic/Fixed Time data**

Substrate depletion: 0.5

Fit limit: 0.95

**Calibrations**

No. of Standards: 0

**Reference Ranges**

Sample type: Serum

Patient...	Min	Max
Female	0.00	0.00
Male	0.00	0.00
Paediatric	0.00	0.00

**Instrument Factor (Y = aX + b)**

a: 1.000 b: 0.000

**Controls**

C1  C2  C3

Dilutions by sample

Formula Builder Print Delete Save

Name	N. of Standards	N. of rep.	Units	Stability
LDH BR	0	1	U/L	

Dil.ratio	Std value	O.D.	Reagent ...

Factor	Factor Min	Factor Max
-6475.000		

KR10282-1/0907