# LDL-CHOLESTEROL

DIRECT

*Enzymatic colorimetric method*

ENDPOINT

<b>REF KR10290</b> 5 x 40 mL <b>CONTENTS</b> R1.Reagent 5 x 30 mL R2.Reagent 5 x 10 mL For <i>in vitro</i> diagnostic use only
---

## SUMMARY

This direct method for quantifying cholesterol in low-density lipoproteins (LDL) is an homogeneous enzymatic test in which the differential precipitation and further sedimentation of the rest of lipoproteins and quilomicrons is avoided. The procedure comprises two steps. In the first step cholesterol in lipoproteins other than LDL in the test sample are decomposed by the simultaneous action of cholesterol esterase (CE) and cholesterol oxidase (CO) at pH 7.0, giving as end products cholestenone and hydrogen peroxide, the latter being decomposed to water and oxygen by catalase. In the second step a surfactant which specifically acts on LDL is added to the reaction product of the first step being the remaining cholesterol quantified by a Trinder's type reaction in which the aniline derivate, HDAOS\*, and 4-aminoantipyrine (4-AA) as a coupling reagent are condensed by the H<sub>2</sub>O<sub>2</sub> in presence of peroxidase (POD) to form a red quinoneimine dye proportional to the concentration of LDL-cholesterol present in the sample.

\* *N*-(2-hydroxy-3-sulfopropyl)-3,5-dimethoxyaniline

## REAGENTS

**R1. Enzyme reagent.** GOOD's buffer 100 mmol/L pH 7.0, CE 5 U/mL, CO 1 U/mL, catalase 200 KU/L, HDAOS 1 mmol/L.

**R2. POD/4-AA reagent.** POD 15 U/mL, 4-AA 4 mmol/L.

## PREPARATION

The Reagents are ready-to-use.

## STORAGE AND STABILITY

 Store at 2-8°C. The Reagents are stable until the expiry date stated on the label. After daily use stored tightly closed and protected from light. On board the reagents are stable 30 days.

### Discard If appear signs of deterioration:

- Presence of particles and turbidity.
- R1.** Blank absorbance (A) at 600 nm > 0.150 in 1cm cuvette.
- R2.** Blank absorbance (A) at 600 nm > 0.100 in 1cm cuvette.

## SAMPLE COLLECTION<sup>1</sup>

Serum, EDTA or heparinized plasma obtained by the patient after an overnight fast. Remove from cells within 3 hours of venipuncture. Samples may be kept at 4-8°C for 2 weeks or at -20°C for 3 months.

## INTERFERENCES

- Lipemia (intralipid >1 g/L) may affect the results.
- Bilirubin (40 mg/dL) does not interfere.
- Hemoglobin (12 g/L) does not interfere.
- Other drugs and substances may interfere<sup>3</sup>.

## INSTRUMENTATION AND MATERIALS

- KROMA analyzer.
- Laboratory equipment.
- iso-Clean Solution. Ref. CT18002.
- Multicalibrator CC/H 10x3 mL Ref. CT19750

## AUTOMATED PROCEDURE

A graphic display pictures the specific sets corresponding to the technical application outlined for this test. Any new application, to the instrument should be validated to confirm that results meet the analytical performance of the method. It is recommended to validate periodically the instrument.

## CALIBRATION

Recalibrate weekly, when a new lot of reagent is used, when control recovery falls out of the expected range or when adjustments are made to the instrument. Two point calibration is recommended (S1: NaCl 9 g/L and S2: Calibrator). A reagent blank should be run daily before sample analysis.

## RESULTS

If results are to be expressed as SI units apply: mg/dL x 0.0259 = mmol/L

## EXPECTED VALUES<sup>2</sup>

Risk group classification according with LDL-C levels.

LDL-Cholesterol levels	RISK
< 100 mg/dL (2.59 mmol/L)	Optimal
100-129 mg/dL (2.59-3.34 mmol/L)	Near optimal
130-159 mg/dL (3.37-4.12 mmol/L).	Borderline
160-189 mg/dL (4.14-4.89 mmol/L)	High
≥ 190 mg/dL (4.92 mmol/L)	Very high

It is recommended that each laboratory establishes its own reference range.

## QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) with assayed values handled as unknowns.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

## DIAGNOSTIC CHARACTERISTICS

Laboratory investigations, epidemiology and genetic factors of hypercholesterolemia indicate that elevated LDH cholesterol is a major cause of coronary heart disease (CHD). The relationship between LDL-cholesterol levels and CHD risk is continuous over a broad range of LDL levels from low to high.

Recent clinical trials robustly show that LDL-lowering therapy reduces risk for CHD. For these reasons the Adult Treatment Panel III continues to identify elevated LDL cholesterol as the primary target of cholesterol-lowering therapy.

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

## PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

## BIBLIOGRAPHY

1. NCCLS Document, "Procedures for the collection of arterial blood specimens". Approved Standard, 3<sup>rd</sup> Ed. (1999).
2. SPECIAL REPORT. Executive Summary of the Third Report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). JAMA. 285 : 2486 (2001).
3. Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AACC Press, 2000.

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