# POTASSIUM BR

*Enzymatic method*

KINETIC

<b>REF KR10330</b> 1 x 40 mL <b>CONTENTS</b> R1. 1 x 32 mL R2. 1 x 8 mL For <i>in vitro</i> diagnostic use only
---

## SUMMARY

Potassium is determined spectrophotometrically through a kinetic coupling assay system using potassium dependent pyruvate kinase. Pyruvate generated is converted to lactate accompanying conversion of NADH analog to NAD analog.

The corresponding decrease of optical density at 380 nm is proportional to the potassium concentration in the serum.

## REAGENT

**R1.** LDH, substrate, NADH analog, azide and stabilizers.

**R2.** Pyruvate Kinase, azide and stabilizers.

## PREPARATION

Ready-to-use.

## STORAGE AND STABILITY

Store at 2-8°C. The Reagents are stable until the expiry date stated on the label. After daily use stored tightly closed and protected from light. On board is stable 30 days. Do not mix reagents of different lots.

### Discard If appear signs of deterioration:

Presence of particles and turbidity.

## SAMPLE COLLECTION

Serum and lithium heparin plasma is recommended. Stable for at least 5 days at 2-8°C. Blood specimens should also be separated from the red cells shortly after collection to prevent any leakage of potassium from the intracellular into the extracellular. Specimens for serum potassium analysis should be free from hemolysis since the high concentration of potassium released from red cells significantly increase the serum levels and this invalidates the test results. Plasma from anticoagulants not containing potassium is also suitable.

## INTERFERENCES

- Na<sup>+</sup> (150 mM), NH<sub>4</sub><sup>+</sup> (0.5 mM), Ca<sup>2+</sup> (7.5 mM), P<sub>i</sub> (2.0 mM), ascorbic acid (10.0 mM), (Zn<sup>2+</sup>) 0.5 mM, Fe<sup>3+</sup> (0.5 mM), Cu<sup>2+</sup> (0.5 mM) does not affect the results.
- Triglycerides (>1000 mg/dL) does not affect the results.
- Hemoglobin (>500 mg/dL) does not affect the results.
- Bilirubin (>20 mg/dL) does not affect the results.
- Other drugs and substances may interfere.

## INSTRUMENTATION AND MATERIALS

- KROMA analyzer.
- Laboratory equipment.
- Cleaning solutions Ref. KR18100, KR18200, KR18300, KR18400.
- Multicalibrator CC/H 5x5 mL Ref. 1975005.

## AUTOMATED PROCEDURE

A graphic display pictures the specific sets corresponding to the technical application outlined for this test.

Any new application, to the instrument should be validated to confirm that results meet the analytical performance of the method.

It is recommended to validate periodically the instrument.

## CALIBRATION

Recalibrate weekly, when a new lot of reagent is used, when control recovery falls out of the expected range or when adjustments are made to the instrument. Two point calibration is recommended (S1: NaCl 9 g/L and S2: Calibrator). A reagent blank should be run daily before sample analysis.

## RESULTS

Samples with concentrations higher than 8.0 mmol/L should be diluted 1:1 with saline and assayed again. Multiply the results by 2.

## EXPECTED VALUES

Adults	3.5 to 5.1 mM (13.7-19.9 mg/dL)
--------	---------------------------------

It is recommended that each laboratory establishes its own reference range.

## QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) Ref. CT19800-CT19850 with assayed values handled as unknowns.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

## DIAGNOSTIC CHARACTERISTICS

Potassium is the principle cation of the intracellular fluid. It is also an important constituent of the extracellular fluid due to its influence on muscle activity. Its intracellular function parallels that of its extracellular function, namely influencing acid-base balance and osmotic pressure, including water retention.

Elevated potassium levels, hyperkalemia, are often associated with renal failure, dehydration shock or adrenal insufficiency. Decreased potassium levels, hypokalemia, are associated with malnutrition, negative nitrogen balance, gastrointestinal fluid losses and hyperactivity of the adrenal cortex.

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

## PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

## BIBLIOGRAPHY

1. Wu, A.H.B., ed. Tietz clinical guide to laboratory tests, 4th edition, p. 880. W.B. Saunders Company, St. Louis (2006).
2. Bergmeyer, H.U., Gawehn, K., and Grassl, M. (1974) in *Methods of Enzymatic Analysis*. Second Edition, Volume I, 509-510, Academic Press, Inc., New York.
3. M.N. Berry, R. D. Mazzachi, M. Pejakovic, and M. J. Peake Enzymatic Determination of Potassium in Serum. *CLIN. CHEM.* 35/5, 817-820 (1989).

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