# SODIUM BR


Enzymatic method

KINETIC

<b>REF KR10340</b> 1 x 40 mL <b>CONTENTS</b> R1. 1 x 26 mL R2. 1 x 13 mL For <i>in vitro</i> diagnostic use only
--

## SUMMARY

Sodium is determined enzymatically via sodium-dependent  $\beta$ -galactosidase activity with ONPG as substrate. The absorbance at 405 nm of the product O-nitrophenyl is proportional to the sodium concentration.

ONPG = *o*-nitrofenil - $\beta$ -D-galactopiranososa

## REAGENT

- R1.** Good's buffer (pH 8.5), Cryptand (>0.4 mM),  $\beta$ -D-galactosidase (<8 U/mL), Proclin 300 (0.02%).
- R2.** Good's buffer (pH 6.5), O-Nitrophenyl  $\beta$ -D- glycoside (>0.5 mM), Proclin 300 (0.02%).

## PREPARATION

Ready-to-use.

## STORAGE AND STABILITY

 Store at 2-8°C. The Reagents are stable until the expiry date stated on the label. After daily use stored tightly closed and protected from light. On board is stable 30 days. Do not mix reagents of different lots.

### Discard if appear signs of deterioration:

Presence of particles and turbidity.

## SAMPLE COLLECTION

Serum. Sodium is stable for at least 24 hours at room temperature and 2 weeks when refrigerated.

## INTERFERENCES

- Triglycerides (>1000 mg/dL) does not affect the results.
- Hemoglobin (>500 mg/dL) does not affect the results.
- Bilirubin (>40 mg/dL) does not affect the results.
- Other drugs and substances may interfere.

## INSTRUMENTATION AND MATERIALS

- KROMA analyzer.
- Laboratory equipment.
- Cleaning solutions Ref. KR18100, KR18200, KR18300, KR18400.
- Multicalibrator CC/H 5x5 mL Ref. 1975005

## AUTOMATED PROCEDURE

A graphic display pictures the specific sets corresponding to the technical application outlined for this test.

Any new application, to the instrument should be validated to confirm that results meet the analytical performance of the method.

It is recommended to validate periodically the instrument.

## CALIBRATION

Recalibrate weekly, when a new lot of reagent is used, when control recovery falls out of the expected range or when adjustments are made to the instrument. Two point calibration is recommended (S1: distilled water and S2: Calibrator). A reagent blank should be run daily before sample analysis.

## RESULTS

Samples with concentrations higher than 180 mmol/L should be diluted 1:1 with saline and assayed again. Multiply the results by 2.

## EXPECTED VALUES

Serum

Adults	135-155 mmol/L (mEq/L)
--------	------------------------

It is recommended that each laboratory establishes its own reference range.

## QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) Ref. CT19800-CT19850 with assayed values handled as unknowns.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

## DIAGNOSTIC CHARACTERISTICS

Sodium is the major cation of extracellular fluid. It plays a central role in the maintenance of the normal distribution of water and the osmotic pressure in the various fluid compartments. The main source of body sodium is sodium chloride contained in ingested foods. Only about one-third of the total body's sodium is contained in the skeleton since most of it is contained in the extracellular body fluids.<sup>1,2</sup>

*Hyponatremia* (low serum sodium level) is found in a variety of conditions including the following: severe polyuria, metabolic acidosis, Addison's disease, diarrhea, and renal tubular disease.

*Hypernatremia* (increased serum sodium level) is found in the following conditions: hyperadrenalism, severe dehydration, diabetic coma after therapy with insulin, excess treatment with sodium salts.

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

## PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

## NOTE

When sodium and potassium are requested together, sodium should be assayed immediately before potassium.

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*In vitro diagnostic reagent for the determination of sodium in serum.*

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**DETERMINACIONES**  
*NUMBER OF TESTS*

**130 tests/kit Code KR10340**

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*(dead volume is not taken in consideration).*

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*Without automatic post-dilution: Up to 180 mmol/L*

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*For additional information read reagent packaging insert.*

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*The reported application is given only as a guideline. Any application to an instrument should be validated to demonstrate that results meet the performance characteristics of the method.*

**General Info**

Name: Sodium Code: Na Barcode: Unit: mmol/l Decimal digits: 0

Visible:  Type: Fixed Time Na+ No. of rgt: 2 Multiply diluted result:

Notes:

---

**Filters**

F1: 405nm F2: not used

**Volumes [microlitres]**

Sample: 8

R1: 200 R2: 100 R3: 0

**Bottle sizes**

R1: 30 ml R2: 20 ml R3: 20 ml

**Incubation/Reading time [sec]**

Substrate/Sample Start

Sample Start: R1,5 -> R2: 360 R1,R2,5 -> R3: 36 Final incub.: 420

Fixed Time second read: 120

**Reference Ranges**

Sample type: Serum

Patient...	Min	Max
Female	0.00	0.00
Male	0.00	0.00
Paediatric	0.00	0.00

**Instrument Factor (Y = aX + b)**

a: 1.000 b: 0.000

**Controls**

C1  C2  C3

**Reagents**

Include blank in calc.

Blank Abs (min; max): 0 0

Reagent linearity: 180

Detection limit: 80

**Calibrations**

No. of Standards: 2

Piecewise Linear Calib. Curve

**Kinetic/Fixed Time data**

Substrate depletion: 0.4

Formula Builder Print Delete Save

Name	N. of Standards	N. of rep.	Units	Stability
Sodium	2	1	mmol/L	

Dil.ratio	Std value	O.D.	Reagent ...
	*		

Factor	Factor Min	Factor Max

\* entered by the user

KR10340-1/0907