PROTEIN

TOTAL


Colorimetric method

ENDPOINT

REF KR10352 8 x 30 mL	REF KR10355 18 x 30 mL
CONTENTS R1.Reagent 8 x 30 mL	CONTENTS R1.Reagent 18 x 30 mL
For <i>in vitro</i> diagnostic use only	

SUMMARY

In the biuret reaction, a chelate is formed between the Cu^{2+} ion and the peptide bonds of the proteins in alkaline solutions to form a violet colored complex whose absorbance is measured photometrically. The intensity of the color produced is proportional to the concentration of protein in the sample.^{1,2}

REAGENT

R1. Biuret reagent. Cupric sulfate 6 mmol/L, sodium-potassium--tartrate 21 mmol/L, potassium iodide 6 mmol/L, sodium hydroxide 0.75 mol/L.
C R:34 S:26-45.

PREPARATION

The Reagent is ready-to-use.

STORAGE AND STABILITY

Store at 2-8°C. The Reagent is stable until the expiry date stated on the label. After daily use stored tightly closed and protected from light. On board this reagent is stable 30 days. Desechar el reactivo cuando presente una absorbancia superior a 0.150 a 540 nm frente agua destilada.

SAMPLE COLLECTION

Serum or heparinized plasma.
Total protein is stable in serum and plasma for 1 week at room temperature, for at least 1 month refrigerated at 2-8°C, and for up to 2 months at -20°C.

INTERFERENCES³

- Lipemia (intralipid) may affect the results.
- Bilirubin (20 mg/dL) does not interfere.
- Hemoglobin may affect the results.
- Other drugs and substances may interfere³.
- Dextrans used as plasma volume expanders for the treatment of low blood pressure, complex with copper and tartrate forming a precipitate.

INSTRUMENTATION AND MATERIALS

- KROMA analyzer.
- Laboratory equipment.
- iso-Clean Solution. Ref. CT18002.
- Multicalibrator CC/H 10x3 mL Ref. CT19750

AUTOMATED PROCEDURE

A graphic display pictures the specific sets corresponding to the technical application outlined for this test.
Any new application, to the instrument should be validated to confirm that results meet the analytical performance of the method.
It is recommended to validate periodically the instrument.

CALIBRATION

Recalibrate weekly, when a new lot of reagent is used, when control recovery falls out of the expected range or when adjustments are made to the instrument. Two point calibration is recommended (S1: NaCl 9 g/L and S2: Calibrator). A reagent blank should be run daily before sample analysis.

RESULTS

Samples with concentrations higher than 12 g/dL should be diluted 1:2 with saline and assayed again. Multiply the results by 2.
If results are to be expressed as SI units apply: g/dL x 10 = g/L.

EXPECTED VALUES⁴

Serum, plasma

Adults	6.6-8.7 g/dL (66-87 g/L)
Prematures	3.6-6.0 g/dL (36-60 g/L)
Newborns	5.3-8.9 g/dL (53-89 g/L)
Pregnancy	Concentration lowers from 69 to 61 g/L

Total serum protein is lower by 4 to 8 g/L with the subject supine than with the subject ambulatory.

Plasma protein is 2 to 4 g/L higher due to the presence of fibrinogen in the sample. It is recommended that each laboratory establishes its own reference range.

QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) with assayed values handled as unknowns. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

DIAGNOSTIC CHARACTERISTICS

The serum content of the soluble proteins, those circulating in extracellular and intracellular fluids, has been used as a marker to aid in clinical diagnosis. The main diagnostic tests are those measuring serum total protein and serum albumin.

Collectively, serum total protein including albumin is mainly involved in the maintenance of normal water distribution between tissues and the blood and responsible for maintaining the oncotic pressure of plasma and is used to transport many substances including macromolecules.

Hiperproteinemia o hiperalbuminemia, usually occurs during multiple mieloma caused by high levels of the monoclonal immunoglobulins, dehydration, excessive water loss, as in severe vomiting, diarrhea, Addison's disease or diabetic acidosis. The *hemoconcentration*, decrease in the volume of plasma water, is reflected as a relative hyperproteinemia since concentration of all the individual plasma proteins are increased to the same degree.

Hypoproteinemia o hypoalbuminemia usually occurs in edema, malnutrition, nephrotic syndrome, malabsorption and severe liver cirrhosis. Since albumin is present in such high concentration low levels of this protein alone may also cause hypoproteinemia.

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

BIBLIOGRAPHY

1. Gornall, A.G., Bardawill, C.S., and David, M.M. J. Biol. Chem. 177 : 751 (1949).
2. Falkner, W.R., and Meites, S. Selected Methods of Clinical Chemistry, 9, 319, AACC., Washington, D.C. (1982).
3. Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AACC Press, 2000.
4. Tietz. N.W. Fundamentals of Clinical Chemistry, p. 940. W.B. Saunders Co. Philadelphia, PA. (1987)

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REFERENCIAS

1. Gornall, A.G., Bardawill, C.S., y David, M.M. J. Biol. Chem. 177 : 751 (1949).
2. Falkner, W.R., and Meites, S. Selected Methods of Clinical Chemistry, 9, 319, AACC., Washington, D.C. (1982).
3. Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AACC Press, 2000.
4. Tietz. N.W. Fundamentals of Clinical Chemistry, p. 940. W.B. Saunders Co. Philadelphia, PA. (1987).

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TOTAL PROTEIN

One reagent method

Total Protein (8 x 30 mL), Code KR10352
Total Protein (18 x 30 mL), Code KR10355

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General Info

Name: Total Protein Code: TP Barcode: c5 Unit: g/dl Decimal digits: 2

Visible: Type: Endpoint Na+: No. of rgt: 1 Multiply diluted result:

Notes:

Filters

F1: 546nm F2: not used

Volumes [microlitres]

Sample: 2

R1: 200 R1 Bottle size: 30 ml

R2: 0 R2 Bottle size: 20 ml

R3: 0 R3 Bottle size: 20 ml

Reagents

Include blank in calc.

Blank Abs (min; max): 0.001 0.8

Reagent linearity: 12

Detection limit: 0.5

Incubation/Reading time [sec]

Substrate/Sample Start: Substrate Start

R1,S -> R2: 300 R1,R2,S -> R3: 36 Final incub.: 300

Reference Ranges

Sample type: Serum

Patient...	Min	Max
Female	0.00	0.00
Male	0.00	0.00
Paediatric	0.00	0.00

Instrument Factor (Y = aX + b)

a: 1.000 b: 0.000

Controls

C1 C2 C3

Calibrations

No. of Standards: 2

Piecewise Linear Calib. Curve

Formula Builder Print Delete Save

Name	N. of Standards	N. of rep.	Units	Stability
Total Protein	2	1	g/dL	

Dil.ratio	Std value	O.D.	Reagent ...
	*		

Factor	Factor Min	Factor Max

* entered by the user

KR10352-55-1/0907