UREA/BUN BR

Urease/GIDH


UV enzymatic method

KINETIC

REF KR10382 5 x 40 mL CONTENTS R1.Reagent 5 x 32 mL R2.Reagent 5 x 8 mL	REF KR10385 11 x 40 mL CONTENTS R1.Reagent 11 x 32 mL R2.Reagent 11 x 8 mL
For <i>in vitro</i> diagnostic use only	

SUMMARY

Urea is hydrolyzed by urease to ammonia and carbon dioxide. The ammonia is converted to glutamate by glutamate dehydrogenase (GIDH) in the presence of NADH and oxoglutarate.^{1,2} The reaction is monitored kinetically at 340 nm by the rate of decrease in absorbance resulting from the oxidation of NADH to NAD⁺, proportional to the concentration of urea present in the sample.

REAGENTS

R1. Buffered Urease/GIDH. TRIS buffer 125 mmol/L pH 7.4, 2-oxoglutarate 10 mmol/L, urease > 140 U/mL, glutamate dehydrogenase > 120 U/mL, Biocides.

R2. Coenzyme. NADH 1.50 mmol/L.

PREPARACIÓN

El Reactivo está listo para su uso.

STORAGE AND STABILITY

Store at 2-8°C. The reagents are stable until the expiry date stated on the label. After daily use stored tightly closed and protected from light. Discard the reagent if presents an absorbance below 1.100 at 340 nm against distilled water. On board the reagents are stable 30 days.

SAMPLE COLLECTION

Serum or heparinized plasma free of hemolysis and urine (see Notes). Other anticoagulants (ammonium heparinate or double oxalate of potassium and ammonium) must not be used.

Urea in serum, plasma or urine is stable 7 days at 2-8°C.

Freeze for longer storage.

INTERFERENCES

- Lipemia (intralipid < 5 g/L) does not interfere.
- Bilirubin (40 mg/dL), hemoglobin (< 4 g/L), do not interfere.
- Other drugs and substances may interfere³.
- Contamination of glassware and water by ammonia, will give falsely elevated results.
- Fluorides used commonly as anticoagulants inhibit the urease of the substrate⁴.

INSTRUMENTATION AND MATERIALS

- KROMA analyzer.
- Laboratory equipment.
- iso-Clean Solution. Ref. CT18002.
- Multicalibrator CC/H 10x3 mL Ref. CT19750

AUTOMATED PROCEDURE

A graphic display pictures the specific sets corresponding to the technical application outlined for this test. Any new application, to the instrument should be validated to confirm that results meet the analytical performance of the method. It is recommended to validate periodically the instrument.

CALIBRATION

Recalibrate weekly, when a new lot of reagent is used, when control recovery falls out of the expected range or when adjustments are made to the instrument. Two point calibration is recommended (S1: NaCl 9 g/L and S2: Calibrator). A reagent blank should be run daily before sample analysis.

RESULTS

Samples with concentrations higher than 500 mg/dL should be diluted 1:2 with saline and assayed again. Multiply the results by 2.

Urine. Dilute the sample 1:50 with distilled water and multiply the result by 50. If results are to be expressed as SI units apply: mg/dL x 0.1665 = mmol/L. To convert urea mass units to those of urea nitrogen apply: mg/dL x 0.467 = mg/dL BUN.

EXPECTED VALUES⁵

Serum, plasma

Newborns (< 10 days)	6.4-53.5 mg/dL (1.1-9.0 mmol/L)
Adults (12-60 years)	15-40 mg/dL (2.5-6.6 mmol/L)

In adults over 60 years of age, the reference interval is 17-50 mg/dL (2.8-8.3 mmol/L) and the concentrations tend to be slightly higher in males than in females.

Urine Adults (normal diet)	26-43 g/24-h (428-714 mmol/24-h)
----------------------------	----------------------------------

A high-protein diet causes significant increases in plasma urea concentrations and urinary excretion.

It is recommended that each laboratory establishes its own reference range.

QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) with assayed values handled as unknowns. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

DIAGNOSTIC CHARACTERISTICS

Azotemia (an abnormal increase in plasma urea level) is seen mainly in renal disorders, dehydration, increase protein catabolism, high-protein diets, or gastrointestinal hemorrhage. There are two types of azotemia. The first, *prerenal azotemia*, is caused by impaired perfusion of the Kidneys due to decreased cardiac output or for any of the former causes. The second, *postrenal azotemia*, is caused by an obstruction in the urine outflow such as nephrolithiasis, prostatism, and tumors of the genitourinary tract. The clinical significance of the urea level in plasma is usually determined in conjugation with the plasma creatinine level. In prerenal azotemia, an increase in the plasma urea level is usually associated with a normal plasma creatinine level, where as in postrenal azotemia, there is an increase in both the urea and the plasma creatinine levels. A decrease in the urea plasma level may be associated with acute dehydration, malnutrition, and pregnancy.

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

NOTES

- Collect a 24-hour urine specimen into a plastic bottle free of preservatives. Keep the sample refrigerated to minimize urea hydrolysis by microorganisms or other agents.

PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

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In vitro diagnostic reagent for the determination of urea in serum, plasma and urine.

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The reported application is given only as a guideline. Any application to an instrument should be validated to demonstrate that results meet the performance characteristics of the method.

General Info

Name: Urea-Bun Code: U Barcode: c7 Unit: mg/dl Decimal digits: 1

Visible: Type: Fixed Time Na+ No. of rgt: 2 Multiply diluted result:

Notes:

Filters

F1: 340nm F2: not used

Volumes [microlitres]

Sample: 2

R1: 160 R2: 40 R3: 0

Bottle sizes

R1: 30 ml R2: 20 ml R3: 20 ml

Incubation/Reading time [sec]

Substrate/Sample Start: Sample Start

R1,S -> R2: 300 R1,R2,S -> R3: 36 Final incub.: 60

Fixed Time second read: 120

Kinetic/Fixed Time data

Substrate depletion: 0.6

Reference Ranges

Sample type: Serum

Patient...	Min	Max
Female	0.00	0.00
Male	0.00	0.00
Paediatric	0.00	0.00

Instrument Factor (Y = aX + b)

a: 1.000 b: 0.000

Controls

C1: C2: C3:

Reagents

Include blank in calc.

Blank Abs (min; max): 1 2

Reagent linearity: 500

Detection limit: 3

Calibrations

No. of Standards: 2

Piecewise Linear Calib. Curve

Formula Builder Print Delete Save

Name	N. of Standards	N. of rep.	Units	Stability
Urea/BUN BR	2	1	md/dL	

Dil.ratio	Std value	O.D.	Reagent ...
	*		

Factor	Factor Min	Factor Max

* entered by the user

KR10382-85-1/0907