# LIPASE

Colorimetric method

ENDPOINT

|   |
|---|
| <b>REF KR10430</b><br>1 x 36 mL<br><b>CONTENTS</b><br>R1. 1 x 30 mL R2. 1 x 6 mL<br>For <i>in vitro</i> diagnostic use only |
|---|

## SUMMARY

The method for the determination of lipase is based on the cleavage of specific chromogenic lipase substrate 1,2-O-dilauryl-rac-glycero-3-glutaric acid-(6'-methyl-resorufin)-ester emulsified in stabilized micro-particles. In the presence of specific activators of pancreatic lipase as colipase, calcium ions and bile acids, the substrate is converted to 1,2-O-dilauryl-rac-glycerol and glutaric acid-6'-methylresorufinester which decomposes spontaneously to glutaric acid and methylresorufin. The increase of absorbance at 580 nm, due to methylresorufin formation, is proportional to the activity of lipase in the sample.

## REAGENT

**R1.** Tris buffer, pH 8.3 40 mmol/L, Colipase  $\geq$  1 mg/L, Desoxycholate  $\geq$  1.8 mM/L, Taurodesoxycholate  $\geq$  7.0 mM/L.

**R2.** Tartrate buffer, pH 4.0 15 mM/L, Lipase Substrate  $\geq$  0.7 mM/L, Calcium Ions  $\geq$  1 mM/L.

**PREPARATION** Ready-to-use.

## STORAGE AND STABILITY

Store at 2-8°C. The Reagents are stable until the expiry date stated on the label. Stability: 90 days at 2-8 °C after opening. After daily use stored tightly closed and protected from light. On board is stable 30 days.

### Discard If appear signs of deterioration:

**R1.** Presence of particles and turbidity.

**R2.** Is a turbid orange-colored micro-emulsion. Discard if it turns red. In some storage conditions (i.e. storage at a temperature lower than the one indicated) a precipitate may appear in the vial that will not influence the reagent performance; however, it is recommended to resuspend the product with a slight rotation of the vial before carrying out the analysis.

## SAMPLE COLLECTION

Serum, EDTA plasma unhemolyzed. Lipase in serum and plasma is stable for 1 week at 2-8°C, and for up to 4 months at -20°C.

## INTERFERENCES

- Triglycerides (>300 mg/dL) may affect the results.
- Hemoglobin (>150 mg/dL) does not affect the results.
- Bilirubin (>20 mg/dL) does not affect the results.
- Other drugs and substances may interfere.

## INSTRUMENTATION AND MATERIALS

- KROMA analyzer.
- Laboratory equipment.
- Cleaning solutions Ref. KR18100, KR18200, KR18300, KR18400.
- Multicalibrator CC/H 10x3 mL Ref. CT19750.

## AUTOMATED PROCEDURE

A graphic display pictures the specific sets corresponding to the technical application outlined for this test.

Any new application, to the instrument should be validated to confirm that results meet the analytical performance of the method.

It is recommended to validate periodically the instrument.

## CALIBRATION

Recalibrate weekly, when a new lot of reagent is used, when control recovery falls out of the expected range or when adjustments are made to the instrument. Two point calibration is recommended (S1: NaCl 9 g/L and S2: Calibrator). A reagent blank should be run daily before sample analysis.

## RESULTS

Samples with concentrations higher than 250 U/L should be diluted 1:10 with saline and assayed again. Multiply the results by 10.

## EXPECTED VALUES

Serum, plasma

|        |          |
|--------|----------|
| Adults | < 38 U/L |
|--------|----------|

It is recommended that each laboratory establishes its own reference range.

## QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) Ref. CT19800-CT19850 with assayed values handled as unknowns.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

## DIAGNOSTIC CHARACTERISTICS

Lipase enzymes are produced in the pancreas and also secreted in small amounts by the salivary glands as well as by gastric, pulmonary and intestinal mucosa. Determination of lipase is used for diagnosis and treatment of diseases of the pancreas such as acute and chronic pancreatitis and obstruction of the pancreatic duct.

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

## PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

## BIBLIOGRAPHY

1. NCCLS. Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard – Fifth Edition (H3-A5). Wayne, PA: The National Committee for Clinical Laboratory Standards, 2003.
2. Kaplan, L.A., Pesce, A.J.: "Clinical Chemistry", Mosby Ed. (1996).
3. EU-Dir 1999/11 Commission Directive of 8 March 1999 adapting to technical progress the principles of Good Laboratory Practice as specified in Council Directive 87/18/EEC.
4. Jakobs, D.S., Kasten, Jr., B.L., DeMott, W.R., Wolfson, W.L.: "Laboratory Test Handbook", Lexi-Comp and Williams & Wilkins Ed. (2nd Edition - 1990).
5. Bonora, R., De Luca, U., Panteghini, M.: "Measurement of pancreatic lipase activity in serum by a kinetic colorimetric assay utilizing a chromogenic substrate reagent". SIBIOC 8-11 october 1996, Pesaro.
6. Neumann, U. et al.: "New substrates for the optical determination of lipase". EP 207252 (1987).

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