FRUCTOSAMINE

Kinetic colorimetric method


FIXED TIME

REF KR10460 2 x 30 mL
CONTENTS R1.Reagent 2 x 30 mL
For <i>in vitro</i> diagnostic use only

SUMMARY

This assay^{1,2} is based on the reducing ability of the ketoamines formed from the non-enzymatic condensation products of glucose and proteins (fructosamines), to reduce tetrazolium salts to a colored dye in alkaline solution.

The oxidized tetrazolium salt, acting also as the reaction starter, is reduced by the fructosamine forming a purple formazan complex.

The rate of color development is proportional to the concentration of fructosamine in the sample.

REAGENT

Fructosamine reagent. Carbonate buffer 0.2 mol/L pH 10.35, NBT 0.5 mmol/L, tensides 2 g/L (w/v). Stabilizers.

Optional. Fructosamine calibrator. Liquid. Glycated secondary standard. Value stated on the label is lot specific.

PREPARATION

The Reagent is ready-to-use.

STORAGE AND STABILITY

 Store at 2-8°C. The Reagent is stable until the expiry date stated on the label. After daily use stored tightly closed and protected from light. On board this reagent is stable 30 days.

Discard the reagent if the blank presents an absorbance over 0.065 at 530 nm against distilled water.

SAMPLE COLLECTION

Serum. Fructosamine in serum is stable for 7 days at 2-8°C.

INTERFERENCES

- Lipids (10 g/L) does not interfere.
- Bilirubin (20 mg/dL) does not interfere.
- Hemoglobin (10 g/L) does not interfere.
- Other drugs and substances may interfere⁴.

INSTRUMENTATION AND MATERIALS

- KROMA analyzer.
- Laboratory equipment.
- Cleaning solutions Ref. KR18100, KR18200, KR18300, KR18400.
- Fructosamine Calibrator CT19460

AUTOMATED PROCEDURE

A graphic display pictures the specific sets corresponding to the technical application outlined for this test.

Any new application, to the instrument should be validated to confirm that results meet the analytical performance of the method.

It is recommended to validate periodically the instrument.

CALIBRATION

Recalibrate weekly, when a new lot of reagent is used, when control recovery falls out of the expected range or when adjustments are made to the instrument. Two point calibration is recommended (S1: Distilled water and S2: Calibrator). A reagent blank should be run daily before sample analysis.

RESULTS

Samples with concentrations higher than 800 µmol/L should be diluted 1:2 with saline and assayed again. Multiply the results by 2.

EXPECTED VALUES³

Serum,

Adults	205-285 µmol/L (glycated albumin)
--------	-----------------------------------

Concentrations are slightly lower (5%) in child. Fructosamine reference values depends on albumin concentration.

Plasmas give lower results than serum.

It is recommended that each laboratory establishes its own reference range.

QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) with assayed values handled as unknowns. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

DIAGNOSTIC CHARACTERISTICS⁴

The fructosamine is formed by non-enzymatic glycosylation of proteins. Since the fructosamine linkages formed are stable fructosamine content is much less variable than glucose content in blood.

Owing to the abnormally high concentration of sugar in blood the amount of fructosamine in serum is increased in diabetes mellitus thus reflecting the degree of glycemic control attained by the diabetic patient and is specially suitable for determining the success of metabolic control by diet and therapeutical treatment over the preceding 10 to 20 days as compared with glycated haemoglobin values which represent a long term (2 to 3 months) index. Higher values of fructosamine are unfailingly indicative of a metabolic poor control.

Disproteïnemic states may produce erroneous results. Severe hypoalbuminemia usually found in conditions such as nephrotic syndrome, may lower serum fructosamine concentrations and not reflect the true degree of the glycemic state. Provided that there is not a severe hypoproteinemia a serum fructosamine concentration will depict the level of glycemia during short periods.

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

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General info

Name: fructosamina Code: Fru Barcode: Unit: $\mu\text{mol/l}$ Decimal digits: 0

Visible: Type: Fixed Time No. of rgt: 1 Multiply diluted result:

Notes:

Filters
F1: 546nm F2: not used

Volumes [microlitres]
Sample: 10
R1: 200 R2: 0 R3: 0
Bottle sizes: R1: 30 ml R2: 20 ml R3: 20 ml

Incubation/Reading time [sec]
Substrate/Sample Start: Substrate Start
R1,S -> R2: 36 R1,R2,S -> R3: 36 Final incub.: 300
Fixed Time second read: 312

Reference Ranges
Sample type: Serum

Patient...	Min	Max
Female	0.00	0.00
Male	0.00	0.00
Paediatric	0.00	0.00

Kinetic/Fixed Time data
Substrate depletion: 0.4

Instrument Factor (Y = aX + b)
a: 1.000 b: 0.000

Reagents
 Include blank in calc.
Blank Abs (min; max): 0 0
Reagent linearity: 800
Detection limit: 0.14

Calibrations
No. of Standards: 2
Piecwise Linear Calib. Curve

Controls
 C1 C2 C3

Dilutions by sample

Formula Builder Print Delete Save

Name	N. of Standards	N. of rep.	Units	Stability
Fructosamine	2	3	$\mu\text{mol/L}$	

Dil.ratio	Std value	O.D.	Reagent ...
	*		

Factor	Factor Min	Factor Max

* entered by the user

KR10460-1/0907