ASLO- Turbidimetric

Latex Turbidimetry

REF KR31100 1 x 40 mL	REF KR31102 5 x 40 mL
CONTENTS R1. Reagent 1 x 32 mL R2. Reagent 1 x 8 mL	CONTENTS R1. Reagent 5 x 32 mL R2. Reagent 5 x 8 mL
For <i>in vitro</i> diagnostic use only	

SUMMARY

The latex particles coated with streptolysin O (SLO) are agglutinated when they react with samples that contain specific antibodies anti-streptolysin O (ASLO). The latex particles agglutination is proportional to the concentration of the ASLO in the sample and can be measured by turbidimetry¹.

REAGENTS

R1. Diluent. Tris buffer, 20 mmol/L, pH 8.2.

R2. Latex. Latex particles coated with estreptolysin O, pH 10.0.

Precautions: The reagents contain sodium azide 0.95 g/L. Avoid any contact with skin or mucous.

PREPARATION

R1. Ready to use.

R2. Ready to use. Mix gently the vial by inversion before use (Note 3).

Calibration curve: Use **ASLO Calibrator Set** (ref. CT39100). The set contains 4 different levels of anti-streptolysin O (ASLO) concentration and it should be used to prepare the calibration curve. Use ClNa 9 g/L or distilled water as calibrator 0. For additional information, read ASLO Calibrator Set Insert.

STORAGE AND STABILITY

-  The reagents will remain stable until the expiration date printed on the label. For optimal stability store **tightly closed** at 2-8°C. Do not use the reagents after the expiration date.
- Reagent deterioration: Presence of particles and turbidity.

SAMPLE COLLECTION

Fresh serum. Stable 7 days at 2-8°C or 3 months at -20°C. Samples with presence of fibrin should be centrifuged before testing. Hemolyzed or contaminated samples are not suitable for testing.

INTERFERENCES

- Bilirubin (40 mg/dL), hemoglobin (12 g/L), lipemia (10 g/L) and rheumatoid factors (800 IU/mL), do not interfere. Other substances may interfere⁸.

ADDITIONAL EQUIPMENT

- KROMA Analyzer.
- Laboratory equipment.
- Cleaning solutions Ref. KR18100, KR18200, KR18300, KR18400.
- Linear ASLO calibrator Ref. 3931105.

AUTHOMATIC TECHNIQUE

For automatic assays the application is included. Any application to an instrument should be validated to demonstrate that results meet the performance characteristics of the method. It is recommended to validate periodically the instrument.

CALIBRATION

Recalibrate weekly, when a new lot of reagent is used, when control recovery falls out of the expected range or when adjustments are made to the instrument. A reagent blank should be run daily before sample analysis.

EXPECTED VALUES⁶⁻⁷

- Adults: Up to 200 IU/mL
- Children (< 2 years): Up to 150 IU/mL.
- Children (school age): Up to 250 IU/mL.

Each laboratory should establish its own reference range.

QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) Ref. 3915010-3915015 with assayed values handled as unknowns.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

DIAGNOSTIC CHARACTERISTICS³⁻⁵

ASLO is a group of specific antibodies developed against and exoenzyme produced by β -hemolytic Streptococci of groups A, C and G. Measuring the ASLO antibodies are useful for the diagnostic of rheumatoid fever, acute glomerulonephritis, bacterial endocarditis and streptococcal infections. Rheumatic fever is an inflammatory disease affecting connective tissue from several parts of human body as (skin, heart, joints etc.) and acute glomerulonephritis is a renal infection that affects mainly to renal glomerulus.

PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

NOTES

- The linearity limit depends on the sample/reagent ratio, as well as the analyzer used. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
- Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.
- For automatic instruments, avoid the presence of bubbles in the reagents that may interfere with the assay results.

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