# $\beta_2$ -microglobulin Turbidimetric

## Latex Turbidimetry

<b>REF KR31160</b> 1 x 40 mL <b>CONTENTS</b> R1. Reagent 1 x 32 mL R2. Reagent 1 x 8 mL	<b>REF KR31162</b> 5 x 40 mL <b>CONTENTS</b> R1. Reagent 5 x 32 mL R2. Reagent 5 x 8 mL
For <i>in vitro</i> diagnostic use only	

### SUMMARY

The latex particles coated with anti-human  $\beta_2$ -microglobulin are agglutinated when they react with samples that containing  $\beta_2$ -microglobulin ( $\beta_2$ -M). The latex particles agglutination is proportional to the concentration of the  $\beta_2$ -M in the sample and can be measured by turbidimetry<sup>1</sup>.

### REAGENTS

**R1 Diluent.** Tris buffer, 20 mmol/L, pH 8.2.

**R2. Latex.** Latex particles coated with goat IgG anti-human  $\beta_2$ -M, pH 8.2.

**Optional:  $\beta_2$ -M Calibrator Set 2 x 1 mL** (ref. CT39160).

**CAL-S Serum Calibrator.** Human serum.

**CAL-U Urine Calibrator.** Human urine.

For additional information, read  $\beta_2$ -M Calibrator Set.

**Precautions:** The reagents contain sodium azide 0.95 g/L. Avoid any contact with skin or mucous.

### PREPARATION

**R1.** Ready to use.

**R2.** Ready to use. Mix gently the vial before use.

### STORAGE AND STABILITY

-  The reagents will remain stable until the expiration date printed on the label. For optimal stability store **tightly closed** at 2-8°C. Do not use the reagents after the expiration date.
- Reagent deterioration: Presence of particles and turbidity.

### SAMPLE COLLECTION

Fresh serum. Stable 7 days at 2-8°C or 3 months at -20°C. Samples with presence of fibrin should be centrifuged before testing. Hemolyzed or contaminated samples are not suitable for testing.

Fresh urine. Adjust urine pH at 7.6 by the addition of 0.4 mol/L  $K_2HPO_4$ .

Stable 2 days at 2-8°C or 2 months at -20°C.

### INTERFERENCES

- Serum:** Bilirubin (40 mg/dL), hemoglobin (12 g/L), lipemia (> 20 g/L) and rheumatoid factors (300 IU/mL), do not interfere.
- Urine:** Hemoglobin (12 g/L), creatinine (300 mg/L), uric acid (500 mg/L) and urea (100 mg/L), do not interfere.  
Other substances may interfere<sup>9</sup>.

### ADDITIONAL EQUIPMENT

- Analyzer KROMA.
- Laboratory equipment.
- Cleaning solutions Ref. KR18100, KR18200, KR18300, KR18400.
- Linear  $\beta_2$ -M Calibrator Set 2 x 1 mL (ref. CT39160).

### AUTHOMATIC TECHNIQUE

For automatic assays the application method is included. Any application to an instrument should be validated to demonstrate that results meet the performance characteristics of the method. It is recommended to validate periodically the instrument.

### CALIBRATION

Recalibrate weekly, when a new lot of reagent is used, when control recovery falls out of the expected range or when adjustments are made to the instrument. A reagent blank should be run daily before sample analysis.

### EXPECTED VALUES<sup>7</sup>

Serum<sup>8</sup>: 1 – 3 mg/L.

Urine<sup>8</sup>: < 0.3 mg/L.

Each laboratory should establish its own reference range.

### QUALITY CONTROL

Control sera are recommended to monitor the performance of manual and automated assay procedures. It is recommended to use Plasma Protein Control N-I (ref: 3915010) and N-II (ref: 3915015).

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

### DIAGNOSTIC CHARACTERISTICS<sup>4-6</sup>

$\beta_2$ -M is a low molecular weight polipeptide (11,800) that is a constituent of HLA<sup>4</sup>. Because approximately 50% of  $\beta_2$ -M is produced by lymphocytes and is freely filtered through glomerular basement membrane,  $\beta_2$ -M levels rise in serum in the presence of glomerular impairment, or lymphocyte activation. Therefore, elevated levels of  $\beta_2$ -M in serum have been reported in disorders such as chronic inflammatory disorders (rheumatoid arthritis, systemic lupus erythematosus, sarcoidosis), viral infection (Epstein Barr virus, cytomegalovirus, influenza A, HIV), lymphoid malignancies (leukemia, malignant lymphoma, myeloma) and cancer (breast, lung, stomach)<sup>4-6</sup>.

As for  $\beta_2$ -M levels in urine, only a trace concentration is detected under normal conditions because up to 99% of  $\beta_2$ -M passed through glomerular is reabsorbed at renal proximal tubules. Therefore  $\beta_2$ -M levels rise significantly with renal tubular dysfunction, or when serum levels are above. The tubular reabsorption threshold (usually > 4 mg/mL)<sup>7</sup>.

### PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

### NOTES

- This method may be used with different instruments. Any application to an instrument should be validated to demonstrate that results meet the performance characteristics of the method. It is recommended to validate periodically the instrument. Contact to the distributor for any question on the application method.
- The linearity limit depends on the sample/reagent ratio, as well as the analyzer used. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
- Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

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Serum/Suero

**General Info**

Name: b2mgser Code: B2SE Barcode: Unit: mg/l Decimal digits: 1

Visible:  Type: Differential - Smp Bll Na+ No. of rgt: 2 Multiply diluted result:

Notes:

---

**Filters**

F1: 546nm F2: not used

**Volumes [microlitres]**

Sample: 3

R1: 240 R2: 60

**Bottle sizes**

R1: 30 ml R2: 20 ml R3: 20 ml

**Incubation/Reading time [sec]**

Substrate/Sample Start: Substrate Start

First incub. R1,R2,S -> R3 Final incub.

120 36 180

**Reference Ranges**

Sample type: Serum

Patient...	Min	Max
Female	1.00	3.00
Male	1.00	3.00
Paediatric	1.00	3.00

**Instrument Factor (Y = aX + b)**

a: 1.000 b: 0.000

**Controls**

C1 C2 C3

**Reagents**

Include blank in calc.

Blank Abs (min; max): 0 0

Reagent linearity: 18

Detection limit: 0

**Calibrations**

No. of Standards: 2

Piecewise Linear Calib. Curve

**Formula Builder** **Print** **Delete** **Save**

Name	N. of Standards	N. of rep.	Units	Stability
$\beta_2$ -M Turbidimetric	2	3	mg/L	

Dil.ratio	Std value	O.D.	Reagent ...
1:1	*		
1:1	0.000		

Factor	Factor Min	Factor Max

\* CAL-S

Urine/Orina

**General Info**

Name: b2mgsor Code: B2OR Barcode: Unit: mg/l Decimal digits: 2

Visible:  Type: Differential - Smp Bl [Na+] No. of rgt: 2 Multiply diluted result:

Notes:

---

**Filters**

F1: 546nm F2: not used

**Volumes [microlitres]**

Sample: 15

R1: 240 R2: 60

**Bottle sizes**

R1: 30 ml R2: 20 ml R3: 20 ml

**Incubation/Reading time [sec]**

Substrate/Sample Start: Substrate Start

First incub. R1,R2,S -> R3 Final incub.

120 36 180

**Reference Ranges**

Sample type: Serum

Patien...	Min	Max
Female	0.00	0.30
Male	0.00	0.30
Paediatric	0.00	0.30

**Instrument Factor (Y = aX + b)**

a: 1.000 b: 0.000

**Controls**

C1 C2 C3

**Reagents**

Include blank in calc.

Blank Abs (min; max): 0 0

Reagent linearity: 3

Detection limit: 0.01

**Calibrations**

No. of Standards: 2

Piecewise Linear Calib. Curve

Formula Builder Print Delete Save

Name	N. of Standards	N. of rep.	Units	Stability
β <sub>2</sub> -M Turbidimetric	2	3	mg/L	

Dil.ratio	Std value	O.D.	Reagent ...
1:1	*		
1:1	0.000		

Factor	Factor Min	Factor Max

\* CAL-U