CRP - Turbidimetric

Latex Turbidimetry

REF KR31200 1 x 40 mL	REF KR31202 5 x 40 mL
CONTENTS R1.Reagent 1 x 32 mL R2. Reagent 1 x 8 mL	CONTENTS R1.Reagent 5 x 32 mL R2.Reagent 5 x 8 mL
For <i>in vitro</i> diagnostic use only	

SUMMARY

The latex particles coated with anti-CRP are agglutinated when they react with samples that contain C-reactive protein (CRP). The latex particles agglutination is proportional to the concentration of the CRP in the sample and can be measured by turbidimetry.^{6,8}

REAGENTS

R1. Diluent. Tris buffer, 20 mmol/L, pH 8.2.

R2. Latex. Latex particles coated with goat anti-human CRP, pH 7.3.

Precautions: The reagents contain sodium azide 0.95 g/L. Avoid any contact with skin or mucous.

PREPARATION

R1. Ready to use.

R2. Ready to use. Shake gently the vial before use.

Calibration curve: Use **CRP Calibrator Set** (ref: CT39200). The set contains 4 different levels of C-Reactive Protein (CRP) concentration and it should be used to prepare the calibration curve. Use ClNa 9 g/L or distilled water as calibrator 0. For additional information, read CRP Calibrator Set Insert.

STORAGE AND STABILITY

-  The reagents will remain stable until the expiration date printed on the label. For optimal stability store **tightly closed** at 2-8°C. Do not use the reagents after the expiration date.
- Reagent deterioration: Presence of particles and turbidity.

SAMPLE COLLECTION

Fresh serum. Stable 7 days at 2-8°C or 3 months at -20°C. Samples with presence of fibrin should be centrifuged before testing. Hemolyzed or contaminated samples are not suitable for testing.

INTERFERENCES

- Bilirubin (40 mg/dL), hemoglobin (12 g/L), lipemia (10 g/L) and rheumatoid factors (800 IU/mL), do not interfere. Other substances may interfere⁹.

ADDITIONAL EQUIPMENT

- Analyzer KROMA.

AUTHOMATIC TECHNIQUE

For automatic assays the application method is included. Any application to an instrument should be validated to demonstrate that results meet the performance characteristics of the method. It is recommended to validate periodically the instrument.

CALIBRATION

Recalibrate weekly, when a new lot of reagent is used, when control recovery falls out of the expected range or when adjustments are made to the instrument. A reagent blank should be run daily before sample analysis.

EXPECTED VALUES^{3,7,8}

Adults: Up to 5 mg/L.

Each laboratory should establish its own reference range.

QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) with assayed values handled as unknown. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

DIAGNOSTIC CHARACTERISTICS¹⁻⁴

La CRP is an acute-phase protein present in normal serum, which increases significantly after most forms of tissue injuries, bacterial and virus infections, inflammation and malignant neoplasia. During tissue necrosis and inflammation resulting from microbial infections, the CRP concentration can rise up to 300 mg/L in 12-24 hours.

PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

NOTES

- The linearity limit depends on the sample/reagent ratio, as well as the analyzer used. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
- Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

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