# FERRITIN- Turbidimetric

*Latex Turbidimetry*

<b>REF KR31600</b> 1 x 40 mL <b>CONTENTS</b> R1.Reagent 1 x 32 mL R2.Reagent 1 x 8 mL	<b>REF KR31602</b> 5 x 40 mL <b>CONTENTS</b> R1.Reagent 5 x 32 mL R2.Reagent 5 x 8 mL
For <i>in vitro</i> diagnostic use only	

## SUMMARY

The latex particles coated with anti human ferritin are agglutinated when they react with samples that contain ferritin. The latex particles agglutination is proportional to the concentration of the ferritin in the sample and can be measured by turbidimetry.<sup>1,2</sup>

## REAGENTS

**R1. Diluent.** Glycine buffer, 20 mmol/L, pH 8.5.

**R2. Latex.** Latex particles coated with polyclonal anti-human ferritin antibodies, pH 8.2.

**Precautions:** The reagents contain sodium azide 0.95 g/L. Avoid any contact with skin or mucous.

## PREPARATION

**R1.** Ready to use.

**R2.** Ready to use. Shake gently the vial before use.

**Calibration curve:** Use **Ferritin Calibrator Set** (ref: CT39400). The set contains 4 different levels of ferritin concentration and it should be used to prepare the calibration curve. Use C1Na 9 g/L or distilled water as calibrator 0. For additional information, read Ferritin Calibrator Set Insert.

## STORAGE AND STABILITY

-  The reagents will remain stable until the expiration date printed on the label. For optimal stability store **tightly closed** at 2-8°C. Do not use the reagents after the expiration date.
- Reagent deterioration: Presence of particles and turbidity.

## SAMPLE COLLECTION

Fresh serum. Stable 7 days at 2-8°C or 3 months at -20°C.

Samples with presence of fibrin should be centrifuged before testing. Hemolyzed or contaminated samples are not suitable for testing.

## INTERFERENCES

- Bilirubin (20 mg/dL), hemoglobin (10 g/L), and rheumatoid factors (600 IU/mL), do not interfere. Lipemia, interferes. Other substances may interfere<sup>8</sup>.

## ADDITIONAL EQUIPMENT

- Analyzer KROMA.

## AUTHOMATIC TECHNIQUE

For automatic assays the application method is included.

Any application to an instrument should be validated to demonstrate that results meet the performance characteristics of the method.

It is recommended to validate periodically the instrument.

## CALIBRATION

Recalibrate weekly, when a new lot of reagent is used, when control recovery falls out of the expected range or when adjustments are made to the instrument. A reagent blank should be run daily before sample analysis.

## EXPECTED VALUES<sup>3-7</sup>

Children: 7 – 140 µg/L

Men: 20 – 250 µg/L

Women: 20 – 200 µg/L

Each laboratory should establish its own reference range.

## QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) with assayed values handled as unknown. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

## DIAGNOSTIC CHARACTERISTICS<sup>3-7</sup>

Ferritin is the major iron storage compound in the body and is considered one of the most reliable indicators of iron status of patients.

A clinical evaluation of serum ferritin is an index of iron stores.

Whereas low serum concentrations of ferritin are always indicative of an iron deficiency, elevated concentrations can occur for variety of reasons. Thus, although elevated concentrations often indicate an excessive iron intake, they are also caused by liver disease, chronic inflammation and malignancies.

Pregnant women, blood donors, hemodialysis patients, adolescents and children are groups particularly at risk. Plasma ferritin is also increased in patients with hemosiderosis or hemochromatosis.

## PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

## NOTES

- The linearity limit depends on the sample/reagent ratio, as well as the analyzer used. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
- Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

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