Microalbumin - Turbidimetric

Latex Turbidimetry

| | |
|---|---|
| REF KR32000 1 x 40 mL CONTENTS R1.Reagent 1 x 32 mL R2.Reagent 1 x 8 mL | REF KR32002 5 x 40 mL CONTENTS R1.Reagent 5 x 32 mL R2.Reagent 5 x 8 mL |
| For <i>in vitro</i> diagnostic use only | |

SUMMARY

Latex particles coated with specific antibodies anti-human albumin are agglutinated when they react with samples that contain albumin. The latex particles agglutination is proportional to the concentration of the albumin in the sample and can be measured by turbidimetry.¹⁻³

REAGENTS

R1 Diluent. Glycin buffer, 100 mmol/L, pH 10.0.

R2. Latex. Latex particles coated with goat IgG anti-human albumin, pH 8.2.

Precautions: The reagents contain sodium azide 0.95 g/L. Avoid any contact with skin or mucous.

PREPARATION

R1. Ready to use.

R2. Ready to use. Mix gently the vial before use.

Calibration curve: Use **Microalbumin Calibrator Set** (ref: CT39500). The set contains 4 different levels of human albumin concentration and it should be used to prepare the calibration curve. Use ClNa 9 g/L or distilled water as calibrator 0.

For additional information, read Microalbumin Calibrator Set Insert.

STORAGE AND STABILITY

-  The reagents will remain stable until the expiration date printed on the label. For optimal stability store **tightly closed** at 2-8°C. Do not use the reagents after the expiration date.
- Reagent deterioration: Presence of particles and turbidity.

SAMPLE COLLECTION

Fresh urine.

It is recommended to adjust the pH at 7.0 with NaOH 1 mol/L. Stable 7 days at 2-8°C when sodium azide 1 g/L is added to prevent contamination. Urine should be centrifuged before testing.

INTERFERENCES

- Bilirubin (10 mg/dL), hemoglobin (12 g/L), urea (100 mg/L) and creatinine (300 mg/L), do not interfere. Other substances may interfere⁸.

ADDITIONAL EQUIPMENT

- Analyzer KROMA.

AUTHOMATIC TECHNIQUE

For automatic assays the application method is included.

Any application to an instrument should be validated to demonstrate that results meet the performance characteristics of the method.

It is recommended to validate periodically the instrument.

CALIBRATION

Recalibrate weekly, when a new lot of reagent is used, when control recovery falls out of the expected range or when adjustments are made to the instrument. A reagent blank should be run daily before sample analysis.

EXPECTED VALUES⁷

Adults: Up to 15 mg/L.

Each laboratory should establish its own reference range.

QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a control with assayed values handled as unknown.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

DIAGNOSTIC CHARACTERISTICS⁴⁻⁶

Microalbuminuria is an increased urinary albumin excretion (UAE) in the range of 20 to 200 µg/min (or 30-300 mg/24h) as a consequence of changes in glomerular permeability⁷.

Increased UAE precedes and is highly predictive of diabetic nephropathy, end-stage renal disease, and proliferative retinopathy in type 1 diabetes. In patients with type 2 diabetes, increased UAE is an independent predictor of progressive renal disease, atherosclerotic disease, and cardiovascular mortality.

In fact, microalbuminuria may show to be a risk factor of cardiovascular disease among otherwise apparently healthy people.

PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

NOTES

- The linearity limit depends on the sample/reagent ratio, as well as the analyzer used. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
- Do not re-use plastic cuvettes, as they may produce erroneous values. Use a new cuvette for each microalbumin assay.
- Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

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